




SURTEX
INSTRUMENTS

Crafted to Excellence




SU



SUTURE

Surtex offers an entire range of instruments for effortless wound closure ranging from needle holders, micro needle holders, ligature needles and many more. Our complete range of high quality suture instruments at your disposal ensures a perfect suture every time.

Our Suture instruments are produced using latest technology to provide outstanding functionality. The diamond dust  DIN-N-DUST coating or tungsten carbide **ULTRAGrip^X** which we use on our forceps and needle holders not only provides a safe grip on tissues and suture material, but also ensure high wear resistance in instruments jaws.

Surtex premium needle holders feature the  DIN-N-DUST coating where the working end of the jaws are coated with a ultra-fine layer of diamond dust the result is a high wear resistance surface which provides firm reliable grip for holding suture needles.

ULTRAGrip^X technology these jaw inserts are engineered by grinding individual teeth into hard, wear resistant tungsten carbide to produce a precise, sharp well-defined point with superb gripping properties.


- Naht
- Sutura
- Suture
- Sutura



Needle Holders
Nadelhalter
Porta-agujas
Porte-aiguilles
Portaghi


Derf Needle Holder
One Fenestrated Jaw

Catalog #	Size
SU-010-12	12 cm (4 3/4")



Converse Needle Holder
Jaw With Longitudinal Serrations

Catalog #	Size
SU-012-12	13 cm (5")



Collier Needle Holder
One Grooved Jaws

Catalog #	Size
SU-016-12	12.5 cm (5")



Webster Needle Holder
Smooth Jaws

Catalog #	Size
SU-020-13	13 cm (5")



Halsey Needle Holder
Smooth Jaws

Catalog #	Size
SU-024-13	13.5 cm (5 1/4")





Needle Holders
Nadelhalter
Porta-agujas
Porte-aiguilles
Portaghi

SU


SUTURE

Halsey Needle Holder
Serrated Jaws



Catalog #	Size
SU-026-13	13.5 cm (5 1/4")

Kilner Needle Holder



Catalog #	Size
SU-029-13	13.5 cm (5 1/4")

Neivert Needle Holder
Smooth Jaws With Special Groove



Catalog #	Size
SU-032-13	12.5 cm (5")

Doyen Needle Holder



Catalog #	Size
SU-034-13	12.5 cm (5")

Hegar-Baumgartner Needle Holder



Catalog #	Size
SU-038-14	14.5 cm (5 3/4")



Needle Holders
Nadelhalter
Porta-agujas
Porte-aiguilles
Portaghi

Brown Needle Holder
Heavy Pattern

Catalog #	Size
SU-040-13	13 cm (5")
SU-040-18	18 cm (7")



Lawrence Needle Holder

Catalog #	Size
SU-042-15	15 cm (6")



Crile-Murray Needle Holder
Slender Pattern

Catalog #	Size
SU-046-15	15 cm (6")



Crile-Wood Needle Holder

Catalog #	Size
SU-050-15	15 cm (6")
SU-050-18	18 cm (7")
SU-050-20	20 cm (8")



Adson Needle Holder
One Fenestrated Jaw

Catalog #	Size
SU-052-18	18 cm (7")





Needle Holders
Nadelhalter
Porta-agujas
Porte-aiguilles
Portaghi

SU

SUTURE

Metzenbaum Needle Holder
One Grooved and Fenestrated Jaws



Catalog #	Size
SU-054-18	18.5 cm (7 1/4")

M.G. Hospital Needle Holder
Heavy Pattern



Catalog #	Size
SU-056-17	17.5 cm (7")

Mayo-Hegar Needle Holder



Catalog #	Size
SU-060-14	14.5 cm (5 3/4")
SU-060-16	16 cm (6 1/4")
SU-060-18	18.5 cm (7 1/4")
SU-060-20	20.5 cm (8")
SU-060-24	24 cm (9 1/2")
SU-060-26	26.5 cm (10 1/2")
SU-060-30	30.5 cm (12")

Mayo-Hegar (Delicate) Needle Holder
Slender Pattern



Catalog #	Size
SU-062-16	16 cm (6 1/4")
SU-062-18	18.5 cm (7 1/4")
SU-062-20	20.5 cm (8")



Needle Holders
Nadelhalter
Porta-agujas
Porte-aiguilles
Portaghi

Sarot Needle Holder



Catalog #	Size
SU-068-18	18 cm (7")
SU-068-26	25.5 cm (10")

Heaney Needle Holder



Catalog #	Size
SU-070-21	20.5 cm (8")

Jameson Needle Holder




Catalog #	Size
SU-072-23	23.5 cm (9 1/4")

Masson Needle Holder



Catalog #	Size
SU-074-27	26.5 cm (10 1/2")

Masson-Luethy Needle Holder



Catalog #	Size
SU-077-26	25.5 cm (10")

INCHES
MILLIMETERS



Finochietto Needle Holder



Catalog #	Size
SU-079-27	26.5 cm (10 1/2")

Wangensteen Needle Holder




Catalog #	Size
SU-082-27	26.5 cm (10 1/2")

Johnson Needle Holder




Catalog #	Size
SU-085-27	27 cm (10 3/4")

Stratte Needle Holder



Catalog #	Size
SU-087-20	20 cm (8")
SU-087-23	23 cm (9")

Olsen-Hegar Needle Holder




Catalog #	Size
SU-090-14	14 cm (5 1/2")
SU-090-18	18.5 cm (7 1/4")



Needle Holders
Nadelhalter
Porta-agujas
Porte-aiguilles
Portaghi

Gillies Needle Holder



Catalog #	Size
SU-094-16	16 cm (6 1/4")

Mathieu Needle Holder



Catalog #	Size
SU-100-14	14 cm (5 1/2")
SU-100-17	17 cm (6 3/4")
SU-100-20	20 cm (8")

Mathieu (Delicate) Needle Holder



Catalog #	Size
SU-102-14	14 cm (5 1/2")
SU-102-17	16.5 cm (6 1/2")

Lichtenberg Needle Holder



Catalog #	Size
SU-110-20	20 cm (8")
SU-110-25	25 cm (9 3/4")


Zweifel Needle Holder



Catalog #	Size
SU-120-20	20 cm (8")
SU-120-23	23 cm (9")



Hoesel Needle Holder



Catalog #	Size
SU-124-25	25 cm (9 3/4")

Crile Needle Holder



Catalog #	Size
SU-130-15	15 cm (6")

Toennis Needle Holder



Catalog #	Size
SU-134-18	18.5 cm (7 1/4")

Langenbeck Needle Holder



Catalog #	Size
SU-140-16	16 cm (6 1/4")
SU-140-18	18 cm (7")
SU-140-20	20 cm (8")

McPherson Micro Needle Holder
Straight, Very Delicate, Without Lock



Catalog #	Size
SU-200-10	10 cm (4")



Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro porte-aiguilles
Micro Portaghi

McPherson Micro Needle Holder
Curved, Very Delicate, Without Lock

Catalog #	Size
SU-201-10	10 cm (4")




McPherson Micro Needle Holder
Straight, Very Delicate, With Lock

Catalog #	Size
SU-202-10	10 cm (4")



McPherson Micro Needle Holder
Curved, Very Delicate, With Lock

Catalog #	Size
SU-203-10	10 cm (4")



Castroviejo Micro Needle Holder
Straight, Delicate

Catalog #	Size
SU-210-14	14 cm (5 1/2")



Castroviejo Micro Needle Holder
Curved, Delicate

Catalog #	Size
SU-211-14	14 cm (5 1/2")





Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro porte-aiguilles
Micro Portaghi

Castroviejo Micro Needle Holder
Straight, Delicate, With Lock

Catalog #	Size
SU-212-14	14 cm (5 1/2")



Castroviejo Micro Needle Holder
Curved, Delicate, With Lock

Catalog #	Size
SU-213-14	14 cm (5 1/2")



Castroviejo Micro Needle Holder
Straight, Standard, Without Lock

Catalog #	Size
SU-220-14	14 cm (5 1/2")



Castroviejo Micro Needle Holder
Curved, Standard, Without Lock

Catalog #	Size
SU-221-14	14 cm (5 1/2")



Castroviejo Micro Needle Holder
Straight, Standard, With Lock

Catalog #	Size
SU-222-14	14 cm (5 1/2")





Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro porte-aiguilles
Micro Portaghi

Castroviejo Micro Needle Holder
Curved, Standard, With Lock



Catalog #	Size
SU-223-14	14 cm (5 1/2")

Jacobson Micro Needle Holder
Straight, Without Lock



Catalog #	Size
SU-230-18	18.5 cm (7 1/4")

Jacobson Micro Needle Holder
Straight, With Lock



Catalog #	Size
SU-231-18	18.5 cm (7 1/4")

Micro Needle Holder
Straight, With Lock



Catalog #	Size
SU-234-21	21 cm (8 1/4")

Yasargil Micro Needle Holder
Straight



Catalog #	Size
SU-238-16	16 cm (6 1/4")

INCHES
MILLIMETERS



Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro porte-aiguilles
Micro Portaghi

Yasargil Micro Needle Holder
Curved

Catalog #	Size
SU-239-16	16 cm (6 1/4")



Kalt Micro Needle Holder
Straight, With Lock

Catalog #	Size
SU-246-14	13.5 cm (5 1/4")




Arruga Micro Needle Holder
Straight, With Lock

Catalog #	Size
SU-250-14	13.5 cm (5 1/4")



Arruga Micro Needle Holder
Curved, With Lock

Catalog #	Size
SU-251-14	13.5 cm (5 1/4")



Castroviejo Micro Needle Holder
Straight, With Lock

Catalog #	Size
SU-260-13	13 cm (5")





Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro porte-aiguilles
Micro Portaghi

Castroviejo Micro Needle Holder
Curved, With Lock



Catalog #	Size
SU-261-13	13 cm (5")

Arruga Needle Holder
Straight, With Ratchet



Catalog #	Size
SU-266-16	16 cm (6 1/4")

Arruga Needle Holder
Curved, With Ratchet



Catalog #	Size
SU-267-16	16 cm (6 1/4")

Micro Needle Holder
Curved, Bayonet Shaped, Smooth Jaws



Catalog #	Size
SU-271-16	16 cm (6 1/4")

Micro Needle Holder
Straight, Bayonet Shaped, Smooth Jaws



Catalog #	Size
SU-274-23	23 cm (9")



Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro porte-aiguilles
Micro Portaghi

Micro Needle Holder

Curved, Bayonet Shaped, Smooth Jaws



Catalog #	Size
SU-275-23	23 cm (9")

Micro Needle Holder

Curved, Bayonet Shaped, Smooth Jaws



Catalog #	Size
SU-277-19	19 cm (7 1/2")

Yasargil Micro Needle Holder

Straight, Bayonet Shaped, Smooth Jaws



Catalog #	Size
SU-280-20	20 cm (8")
SU-280-23	23 cm (9")

Yasargil Micro Needle Holder

Curved, Bayonet Shaped, Smooth Jaws



Catalog #	Size
SU-281-20	20 cm (8")
SU-281-23	23 cm (9")

Barraquer-Troutman Micro Needle Holder

Straight, Extra Delicate, Round Handle



Catalog #	Size
SU-300-10	10 cm (4")



Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro porte-aiguilles
Micro Portaghi

Barraquer-Troutman Micro Needle Holder
Curved, Extra Delicate, Round Handle



Catalog #	Size
SU-301-10	10.5 cm (4")

Barraquer-Troutman Micro Needle Holder
Straight, Very Delicate, Round Handle, With Lock



Catalog #	Size
SU-310-10	10.5 cm (4")
SU-310-11	11.5 cm (4½")

Barraquer-Troutman Micro Needle Holder
Curved, Very Delicate, Round Handle, With Lock



Catalog #	Size
SU-311-10	10.5 cm (4")
SU-311-11	11.5 cm (4 ½")

Barraquer Micro Needle Holder
Curved, Delicate, Round Handle



Catalog #	Size
SU-313-11	11.5 cm (4 ½")

Barraquer Micro Needle Holder
Curved, Delicate, Round Handle, With Lock



Catalog #	Size
SU-315-11	11.5 cm (4 ½")



Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro porte-aiguilles
Micro Portaghi

Titan Barraquer Micro Needle Holder

Curved, Delicate, Round Handle, With Lock, Titanium

Catalog #	Size
SU-315-11 T	11.5 cm (4 1/2")



Barraquer Micro Needle Holder

Straight, Very Delicate, Round Handle

Catalog #	Size
SU-320-13	13 cm (5")



Barraquer Micro Needle Holder

Curved, Very Delicate, Round Handle

Catalog #	Size
SU-321-13	13 cm (5")



Barraquer Micro Needle Holder

Straight, Very Delicate, Round Handle, With Lock

Catalog #	Size
SU-322-13	13 cm (5")



Barraquer Micro Needle Holder

Curved, Very Delicate, Round Handle, With Lock

Catalog #	Size
SU-323-13	13 cm (5")






Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro porte-aiguilles
Micro Portaghi

Barraquer Micro Needle Holder
Curved, Very Delicate, Round Handle



Catalog #	Size
SU-331-14	14 cm (5 1/2")

Barraquer Micro Needle Holder
Curved, Very Delicate, Round Handle, Titanium




Catalog #	Size
SU-331-14 T	14 cm (5 1/2")

Barraquer Micro Needle Holder
Curved, Very Delicate, Round Handle, With Lock



Catalog #	Size
SU-333-14	14 cm (5 1/2")

Barraquer Micro Needle Holder
Curved, Very Delicate, Round Handle, With Lock, Titanium



Catalog #	Size
SU-333-14 T	14 cm (5 1/2")

Troutman Micro Needle Holder
Strongly Curved, Delicate, Round Handle



Catalog #	Size
SU-341-12	12 cm (4 3/4")



Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro porte-aiguilles
Micro Portaghi

Troutman Micro Needle Holder

Strongly Curved, Delicate, Round Handle, With Lock



Catalog #	Size
SU-343-12	12 cm (4 3/4")

Reill Micro Needle Holder

Straight



Catalog #	Size
SU-350-15	15 cm (6")

Reill Micro Needle Holder

Curved



Catalog #	Size
SU-351-15	15 cm (6")

Micro Needle Holder

Straight, With Lock



Catalog #	Size
SU-352-15	15 cm (6")

Micro Needle Holder

Curved, With Lock




Catalog #	Size
SU-353-15	15 cm (6")




Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro porte-aiguilles
Micro Portaghi

Micro Needle Holder
Straight



Catalog #	Size
SU-354-18	18.5 cm (7 1/4")

Micro Needle Holder
Curved



Catalog #	Size
SU-355-18	18.5 cm (7 1/4")

Micro Needle Holder
Straight, With Lock



Catalog #	Size
SU-356-18	18.5 cm (7 1/4")

Micro Needle Holder
Curved, With Lock



Catalog #	Size
SU-357-18	18.5 cm (7 1/4")

Micro Needle Holder
Straight, With Lock



Catalog #	Size
SU-360-20	20 cm (8")

INCHES
MILLIMETERS



Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro porte-aiguilles
Micro Portaghi

Micro Needle Holder
Straight, Delicate, With Lock



Catalog # Size
SU-362-23 23 cm (9")

Diam-n-Dust™ Castroviejo Micro Needle Holder
Straight, Extra Delicate



Catalog # Size
SU-400-14 14 cm (5 1/2")



Diam-n-Dust™ Castroviejo Micro Needle Holder
Curved, Extra Delicate



Catalog # Size
SU-401-14 14 cm (5 1/2")



Diam-n-Dust™ Castroviejo Micro Needle Holder
Straight, Extra Delicate, With Lock



Catalog # Size
SU-402-14 14 cm (5 1/2")



Diam-n-Dust™ Castroviejo Micro Needle Holder
Curved, Extra Delicate, With Lock



Catalog # Size
SU-403-14 14 cm (5 1/2")





Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro porte-aiguilles
Micro Portaghi

Diam-n-Dust™ Castroviejo Micro Needle Holder
Straight, Very Delicate



DIAM-N-DUST

Catalog #	Size
SU-404-14	14 cm (5 1/2")



Diam-n-Dust™ Castroviejo Micro Needle Holder
Curved, Very Delicate



DIAM-N-DUST

Catalog #	Size
SU-405-14	14 cm (5 1/2")



Diam-n-Dust™ Castroviejo Micro Needle Holder
Straight, Very Delicate, With Lock



DIAM-N-DUST

Catalog #	Size
SU-406-14	14 cm (5 1/2")



Diam-n-Dust™ Castroviejo Micro Needle Holder
Curved, Very Delicate, With Lock



DIAM-N-DUST

Catalog #	Size
SU-407-14	14 cm (5 1/2")



Diam-n-Dust™ Castroviejo Micro Needle Holder
Straight, Delicate



DIAM-N-DUST

Catalog #	Size
SU-410-18	18 cm (7")





Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro porte-aiguilles
Micro Portaghi

Diam-n-Dust™ Castroviejo Micro Needle Holder
Curved, Delicate



Catalog #	Size
SU-411-18	18 cm (7")



Diam-n-Dust™ Castroviejo Micro Needle Holder
Straight, Delicate, With Lock



Catalog #	Size
SU-412-18	18 cm (7")



Diam-n-Dust™ Castroviejo Micro Needle Holder
Curved, Delicate, With Lock



Catalog #	Size
SU-413-18	18 cm (7")



Diam-n-Dust™ Micro Needle Holder
Straight, Without Lock



Catalog #	Size
SU-420-14	14 cm (5 1/2")
SU-420-18	18 cm (7")
SU-420-21	21 cm (8 1/4")
SU-420-23	23 cm (9")



Diam-n-Dust™ Micro Needle Holder
Straight, With Lock



Catalog #	Size
SU-422-14	14 cm (5 1/2")
SU-422-18	18 cm (7")
SU-422-21	21 cm (8 1/4")
SU-422-23	23 cm (9")





Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro porte-aiguilles
Micro Portaghi

Diam-n-Dust™ Micro Needle Holder
Straight, Round Handle



DIAM-N-DUST

Catalog #	Size
SU-430-14	14 cm (5 1/2")



Diam-n-Dust™ Titan Micro Needle Holder
Straight, Round Handle, Titanium



DIAM-N-DUST

Catalog #	Size
SU-430-14 T	14 cm (5 1/2")



Diam-n-Dust™ Micro Needle Holder
Straight, Round Handle



DIAM-N-DUST

Catalog #	Size
SU-430-16	16 cm (6 1/4")



Diam-n-Dust™ Titan Micro Needle Holder
Straight, Round Handle, Titanium



DIAM-N-DUST

Catalog #	Size
SU-430-16 T	16 cm (6 1/4")



Diam-n-Dust™ Micro Needle Holder
Straight, Round Handle



DIAM-N-DUST

Catalog #	Size
SU-430-18	18 cm (7")





Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro porte-aiguilles
Micro Portaghi

Diam-n-Dust™ Titan Micro Needle Holder
Straight, Round Handle, Titanium



DIAM-N-DUST

Catalog # Size
SU-430-18 T 18 cm (7")



Diam-n-Dust™ Micro Needle Holder
Straight, Round Handle



DIAM-N-DUST

Catalog # Size
SU-430-21 21 cm (8 1/4")



Diam-n-Dust™ Titan Micro Needle Holder
Straight, Round Handle, Titanium



DIAM-N-DUST

Catalog # Size
SU-430-21 T 21 cm (8 1/4")



Diam-n-Dust™ Micro Needle Holder
Straight, Round Handle



DIAM-N-DUST

Catalog # Size
SU-430-23 23 cm (9")



Diam-n-Dust™ Micro Needle Holder
Straight, Round Handle, With Lock



DIAM-N-DUST

Catalog # Size
SU-432-14 14 cm (5 1/2")





Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro porte-aiguilles
Micro Portaghi

Diam-n-Dust™ Titan Micro Needle Holder
Straight, Round Handle, With Lock, Titanium



DIAM-N-DUST

Catalog #	Size
SU-432-14 T	14 cm (5 1/2")



Diam-n-Dust™ Micro Needle Holder
Straight, Round Handle, With Lock



DIAM-N-DUST

Catalog #	Size
SU-432-16	16 cm (6 1/4")




Diam-n-Dust™ Titan Micro Needle Holder
Straight, Round Handle, With Lock, Titanium



DIAM-N-DUST

Catalog #	Size
SU-432-16 T	16 cm (6 1/4")



Diam-n-Dust™ Micro Needle Holder
Straight, Round Handle, With Lock



DIAM-N-DUST

Catalog #	Size
SU-432-18	18 cm (7")



Diam-n-Dust™ Titan Micro Needle Holder
Straight, Round Handle, With Lock, Titanium



DIAM-N-DUST

Catalog #	Size
SU-432-18 T	18 cm (7")





Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro porte-aiguilles
Micro Portaghi

Diam-n-Dust™ Micro Needle Holder
Straight, Round Handle, With Lock



DIA.M-N-DUST

Catalog #	Size
SU-432-21	21 cm (8 1/4")



Diam-n-Dust™ Titan Micro Needle Holder
Straight, Round Handle, With Lock, Titanium



DIA.M-N-DUST

Catalog #	Size
SU-432-21 T	21 cm (8 1/4")



Diam-n-Dust™ Micro Needle Holder
Straight, Round Handle, With Lock



DIA.M-N-DUST

Catalog #	Size
SU-432-23	23 cm (9")



Diam-n-Dust™ Micro Needle Holder
Straight, Heavy Pattern, Round Handle



DIA.M-N-DUST

Catalog #	Size
SU-440-19	19 cm (7 1/2")
SU-440-23	23 cm (9")
SU-440-25	25 cm (9 3/4")



Diam-n-Dust™ Micro Needle Holder
Straight, Heavy Pattern, Round Handle, With Lock



DIA.M-N-DUST

Catalog #	Size
SU-442-19	19 cm (7 1/2")
SU-442-23	23 cm (9")
SU-442-25	25 cm (9 3/4")





Needle Holders
Nadelhalter
Porta-agujas
Porte-aiguilles
Portaghi

Diam-n-Dust™ Castroviejo Micro Needle Holder
Straight, With Lock




Catalog #	Size
SU-450-12	12.5 cm (5")

Diam-n-Dust™ Hepp-Scheidel Micro Needle Holder
Straight, With Lock



Catalog #	Size
SU-452-18	18.5 cm (7 1/4")

Diam-n-Dust™ DeBakey Needle Holder



Catalog #	Size
SU-460-15	15 cm (6")
SU-460-18	18 cm (7")
SU-460-20	20 cm (8")
SU-460-23	23 cm (9")

Diam-n-Dust™ Ryder Vascular Needle Holder



Catalog #	Size
SU-470-13	13 cm (5")
SU-470-15	15 cm (6")
SU-470-18	18 cm (7")
SU-470-22	22 cm (8 3/4")
SU-470-25	25 cm (9 3/4")



Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro porte-aiguilles
Micro Portaghi

UltraGrip™ TC Barraquer Micro Needle Holder
Curved, Smooth Jaws, Round Handle



ULTRAGrip X™



Catalog #	Size
SU-500-13	13 cm (5")

UltraGrip™ TC Barraquer Micro Needle Holder
Curved, Smooth Jaws, Round Handle , With Lock



ULTRAGrip X™



Catalog #	Size
SU-502-13	13 cm (5")

UltraGrip™ TC Castroviejo Micro Needle Holder
Straight, With Lock



ULTRAGrip X™



Catalog #	Size
SU-510-14	14.5 cm (5 3/4")

UltraGrip™ TC Castroviejo Micro Needle Holder
Curved, With Lock



ULTRAGrip X™



Catalog #	Size
SU-511-14	14.5 cm (5 3/4")

UltraGrip™ TC Castroviejo Micro Needle Holder
Straight, Smooth Jaws, With Lock



ULTRAGrip X™



Catalog #	Size
SU-512-14	14.5 cm (5 3/4")



Micro Needle Holders
Mikro-Nadelhalter
Micro Porta-agujas
Micro porte-aiguilles
Micro Portaghi

UltraGrip™ TC Castroviejo Micro Needle Holder
Curved, Smooth Jaws, With Lock



ULTRAGrip™ X™

Catalog #	Size
SU-513-14	14.5 cm (5 3/4")



UltraGrip™ TC Micro Needle Holder
Straight, With Lock



ULTRAGrip™ X™

Catalog #	Size
SU-520-18	18 cm (7")



UltraGrip™ TC Micro Needle Holder
Curved, With Lock



ULTRAGrip™ X™

Catalog #	Size
SU-521-18	18 cm (7")



UltraGrip™ TC Jacobson Micro Needle Holder
Straight, With Lock



ULTRAGrip™ X™

Catalog #	Size
SU-530-21	21.5 cm (8 1/2")



UltraGrip™ TC Jacobson Micro Needle Holder
Curved, With Lock



Catalog #	Size
SU-531-21	21.5 cm (8 1/2")





UltraGrip™ TC Castroviejo Micro Needle Holder
Smooth Jaws, With Lock

ULTRAGripX™

Catalog #	Size
SU-540-13	12.5 cm (5")



UltraGrip™ TC Hepp-Scheidel Micro Needle Holder
Smooth Jaws, With Lock

ULTRAGripX™

Catalog #	Size
SU-544-18	18.5 cm (7 1/4")



UltraGrip™ TC Diethrich Micro Needle Holder

ULTRAGripX™


Catalog #	Size
SU-550-16	16 cm (6 1/4")
SU-550-18	18 cm (7")
SU-550-20	20 cm (8")
SU-550-23	23 cm (9")
SU-550-26	26 cm (10 1/2")



UltraGrip™ TC Arruga Needle Holder
Straight

ULTRAGripX™

Catalog #	Size
SU-560-16	16 cm (6 1/4")



UltraGrip™ TC Arruga Needle Holder
Curved

Catalog #	Size
SU-561-16	16 cm (6 1/4")






Needle Holders
Nadelhalter
Porta-agujas
Porte-aiguilles
Portaghi

UltraGrip™ TC Webster Needle Holder
With Automatic Release Ratchet

ULTRAGrip X™


Catalog #	Size
SU-564-13	13 cm (5")



UltraGrip™ TC Crile-Wood Needle Holder
With Automatic Release Ratchet

ULTRAGrip X™

Catalog #	Size
SU-566-15	15 cm (6")



UltraGrip™ TC Webster Needle Holder
Smooth Jaws

ULTRAGrip X™

Catalog #	Size
SU-570-13	13 cm (5")



UltraGrip™ TC Kilner Needle Holder
S Shaped

ULTRAGrip X™

Catalog #	Size
SU-573-13	13 cm (5")



UltraGrip™ TC Derf Needle Holder

ULTRAGrip X™

Catalog #	Size
SU-576-12	12 cm (4 3/4")





UltraGrip™ TC Hegar-Baumgartner Needle Holder

ULTRAGripX™

Catalog #	Size
SU-580-13	13 cm (5")
SU-580-14	14.5 cm (5")



**UltraGrip™ TC Halsey Needle Holder
Smooth Jaw**

ULTRAGripX™

Catalog #	Size
SU-584-13	13 cm (5")



UltraGrip™ TC Halsey Needle Holder

ULTRAGripX™


Catalog #	Size
SU-586-13	13 cm (5")



**UltraGrip™ TC Webster Needle Holder
Smooth Jaw**

ULTRAGripX™


Catalog #	Size
SU-588-12	12 cm (4 3/4")



UltraGrip™ TC Converse Needle Holder

ULTRAGripX™

Catalog #	Size
SU-590-12	12.5 cm (5")





Needle Holders
Nadelhalter
Porta-agujas
Porte-aiguilles
Portaghi

UltraGrip™ TC Crile-Wood Needle Holder

ULTRAGrip™


Catalog #	Size
SU-594-15	15 cm (6")
SU-594-18	18 cm (7")
SU-594-20	20 cm (8")



UltraGrip™ TC Sarot Needle Holder

ULTRAGrip™


Catalog #	Size
SU-598-18	18 cm (7")
SU-598-20	20 cm (8")
SU-598-26	26 cm (10 1/4")



UltraGrip™ TC Mayo-Hegar Needle Holder

ULTRAGrip™

Catalog #	Size
SU-600-14	14.5 cm (5 3/4")
SU-600-16	16 cm (6 1/4")
SU-600-18	18 cm (7")
SU-600-20	20 cm (8")
SU-600-24	24 cm (9 1/2")
SU-600-26	26 cm (10 1/4")
SU-600-30	30 cm (11 3/4")



UltraGrip™ TC Hegar Needle Holder

ULTRAGrip™

Catalog #	Size
SU-602-20	20 cm (8")
SU-602-24	23.5 cm (9 1/4")



INCHES
MILLIMETERS



UltraGrip™ TC Heaney Needle Holder

ULTRAGripX™

Catalog #	Size
SU-605-20	20 cm (8")
SU-605-24	24 cm (9 ½")
SU-605-26	26 cm (10 ¼")



UltraGrip™ TC Bozemann Needle Holder
S Shaped

ULTRAGripX™


Catalog #	Size
SU-611-20	20.5 cm (8")
SU-611-24	24 cm (9 ½")
SU-611-26	26.5 cm (10 ½")
SU-611-28	29 cm (11 ½")



UltraGrip™ TC Stratte Needle Holder

ULTRAGripX™

Catalog #	Size
SU-615-21	21 cm (8 ¼")
SU-615-23	23 cm (9")



UltraGrip™ TC Masson Needle Holder

ULTRAGripX™

Catalog #	Size
SU-620-27	26.5 cm (10 ½")





Needle Holders
Nadelhalter
Porta-agujas
Porte-aiguilles
Portaghi

UltraGrip™ TC Wangensteen Needle Holder

ULTRAGrip X™


Catalog #	Size
SU-624-27	27 cm (10 3/4")



UltraGrip™ TC Needle Holder

ULTRAGrip X™

Catalog #	Size
SU-625-27	27 cm (10 3/4")



UltraGrip™ TC Finochietto Needle Holder

ULTRAGrip X™

Catalog #	Size
SU-627-27	26.5 cm (10 1/2")



UltraGrip™ TC Naunton-Morgan Needle Holder
Compound Action

ULTRAGrip X™


Catalog #	Size
SU-630-26	26 cm (10 1/4")



UltraGrip™ TC DeBakey Needle Holder

ULTRAGrip X™

Catalog #	Size
SU-650-18	18 cm (7")
SU-650-20	20 cm (8")
SU-650-23	23 cm (9")
SU-650-26	26 cm (10 1/4")
SU-650-31	31 cm (12 1/4")





UltraGrip™ TC Ryder Needle Holder

ULTRAGripX™


Catalog #	Size
SU-660-13	12.5 cm (5")
SU-660-15	15 cm (6")
SU-660-18	18 cm (7")
SU-660-20	20 cm (8")
SU-660-23	23 cm (9")
SU-660-26	26 cm (10 1/4")



UltraGrip™ TC Ryder-Vascular Needle Holder

ULTRAGripX™


Catalog #	Size
SU-664-13	12.5 cm (5")
SU-664-15	15 cm (6")
SU-664-18	18 cm (7")
SU-664-20	20 cm (8")
SU-664-23	22.5 cm (8 3/4")
SU-664-26	26 cm (10 1/4")



UltraGrip™ TC Martin Needle Holder

ULTRAGripX™

Catalog #	Size
SU-670-15	15 cm (6")
SU-670-17	17 cm (6 3/4")
SU-670-20	20 cm (8")






Needle Holders
Nadelhalter
Porta-agujas
Porte-aiguilles
Portaghi

UltraGrip™ TC Hegar-Vascular Needle Holder

ULTRAGrip X™

Catalog #	Size
SU-674-14	14 cm (5 1/2")
SU-674-18	18.5 cm (7 1/4")
SU-674-20	20.5 cm (8")
SU-674-23	23 cm (9")



UltraGrip™ TC Micro Vascular Needle Holder

ULTRAGrip X™


Catalog #	Size
SU-678-14	14 cm (5 1/2")
SU-678-18	18 cm (7")
SU-678-20	20 cm (8")
SU-678-23	23 cm (9")



UltraGrip™ TC Martini Needle Holder
Smooth Jaw

ULTRAGrip X™

Catalog #	Size
SU-680-14	14 cm (5 1/2")
SU-680-16	16 cm (6 1/4")
SU-680-18	18 cm (7")
SU-680-20	20 cm (8")



UltraGrip™ TC Martini Needle Holder

ULTRAGrip X™

Catalog #	Size
SU-682-14	14 cm (5 1/2")
SU-682-16	16 cm (6 1/4")
SU-682-18	18 cm (7")
SU-682-20	20 cm (8")





UltraGrip™ TC Olsen-Hegar Needle Holder

ULTRAGrip X™

Catalog #	Size
SU-690-14	14.5 cm (5 ¾")
SU-690-17	17 cm (6 ¾")
SU-690-18	19 cm (7 ½")



UltraGrip™ TC Foster-Gillies Needle Holder
Smooth Jaws

ULTRAGrip X™

Catalog #	Size
SU-693-11	11 cm (4 ¼")



UltraGrip™ TC Gillies-Mini Needle Holder

ULTRAGrip X™


Catalog #	Size
SU-695-15	15 cm (6")



UltraGrip™ TC Gillies Needle Holder

ULTRAGrip X™

Catalog #	Size
SU-697-16	16 cm (6 ¼")



UltraGrip™ TC Mathieu Needle Holder

ULTRAGrip X™

Catalog #	Size
SU-700-14	14 cm (5 ½")
SU-700-17	17 cm (6 ¾")
SU-700-20	20 cm (8")





Needle Holders
Nadelhalter
Porta-agujas
Porte-aiguilles
Portaghi

UltraGrip™ TC Mathieu Needle Holder

ULTRAGrip X™

Catalog #	Size
SU-702-14	14 cm (5 1/2")
SU-702-17	17 cm (6 3/4")
SU-702-20	20 cm (8")



**UltraGrip™ TC Mathieu Needle Holder
With Interior Ratchet**

ULTRAGrip X™

Catalog #	Size
SU-706-17	17 cm (6 3/4")
SU-706-20	20 cm (8")
SU-706-23	23 cm (9")



UltraGrip™ TC Lichtenberg Needle Holder

ULTRAGrip X™

Catalog #	Size
SU-710-20	20 cm (8")



**UltraGrip™ TC Zweifel Needle Holder
With Hinged Ratchet**

ULTRAGrip X™

Catalog #	Size
SU-720-20	20.5 cm (8")
SU-720-23	23 cm (9")





Needle Holder, Ligature Needles
Nadelhalter, Unterbindungsnadeln
Porta-agujas, Agujas para ligaduras
Porte-aiguilles, Aiguilles à ligatures
Portaghi, Aghi per legatura

UltraGrip™ TC Toennis Needle Holder

ULTRAGrip™ X

Catalog #	Size
SU-730-18	18.5 cm (7 1/4")



UltraGrip™ TC Langenbeck Needle Holder

ULTRAGrip™ X

Catalog #	Size
SU-740-18	18 cm (7")
SU-740-20	20.5 cm (8")
SU-740-23	23 cm (9")




Reverdin Ligature Needle Short

Catalog #	Size
SU-801-15	14.5 cm (5 3/4")




Reverdin Ligature Needle Short

Catalog #	Size
SU-803-15	14.5 cm (5 3/4")



Reverdin Ligature Needle Short


Catalog #	Size
SU-805-15	14.5 cm (5 3/4")






Ligature Needles
Unterbindungsnadeln
Agujas para ligaduras
Aiguilles à ligatures
Aghi per legatura

Reverdin Ligature Needle
Short




Catalog #	Size
SU-807-15	14.5 cm (5 3/4")

Reverdin Ligature Needle
Medium




Catalog #	Size
SU-811-19	19.5 cm (7 3/4")

Reverdin Ligature Needle
Medium




Catalog #	Size
SU-813-19	19.5 cm (7 3/4")

Reverdin Ligature Needle
Medium



Catalog #	Size
SU-815-19	19.5 cm (7 3/4")

Reverdin Ligature Needle
Medium




Catalog #	Size
SU-817-19	19.5 cm (7 3/4")

INCHES
MILLIMETERS



Ligature Needles
Unterbindungsnadeln
Agujas para ligaduras
Aiguilles à ligatures
Aghi per legatura


Reverdin Ligature Needle
Long



Catalog # Size
SU-821-23 23 cm (9")




Reverdin Ligature Needle
Long



Catalog # Size
SU-823-23 23 cm (9")




Reverdin Ligature Needle
Long



Catalog # Size
SU-825-23 23 cm (9")



Reverdin Ligature Needle
Long



Catalog # Size
SU-827-23 23 cm (9")



Young Ligature Needle Holder Complete
With Needle Ref: SU-831-01




Catalog # Size
SU-831-00 24.5 cm (9 3/4")



Ligature Needles
Unterbindungsnadeln
Agujas para ligaduras
Aiguilles à ligatures
Aghi per legatura


Young Ligature Needle Only
Fig. 1

Catalog #	Size
SU-831-01	N/A




Young Ligature Needle Only
Fig. 2

Catalog #	Size
SU-831-02	N/A




Young Ligature Needle Only
Fig. 3

Catalog #	Size
SU-831-03	N/A



Schmieden-Dick Ligature Needle
Malleable

Catalog #	Size
SU-835-32	29 cm (11 1/2")



Schmieden Ligature Needle

Catalog #	Size
SU-837-21	24 cm (9 1/2")





Ligature Needles
Unterbindungsnadeln
Agujas para ligaduras
Aiguilles à ligatures
Aghi per legatura

Ligature Needle
Malleable



Catalog #	Size
SU-841-20	20.5 cm (8")

Jonesco Ligature Needle




Catalog #	Size
SU-843-18	18.5 cm (7 1/4")

Claus Van Eicken Ligature Needle




Catalog #	Size
SU-845-25	25 cm (9 3/4")

Ligature Needle



Catalog #	Size
SU-849-15	15 cm (6")

Syme Ligature Needle




Catalog #	Size
SU-851-17	17 cm (6 3/4")



Ligature Needles
Unterbindungsnadeln
Agujas para ligaduras
Aiguilles à ligatures
Aghi per legatura


Cooper Ligature Needle



Catalog #	Size
SU-853-19	18.5 cm (7 1/4")



Kocher Ligature Needle



Catalog #	Size
SU-855-20	20.5 cm (8")



Kronecker Ligature Needle

Delicate For Right Hand



Catalog #	Size
SU-861-13	12.5 cm (5")



Kronecker Ligature Needle

Delicate For Left Hand




Catalog #	Size
SU-863-13	12.5 cm (5")



Deschamps Ligature Needle

Blunt For Right Hand



Catalog #	Size
SU-871-24	24 cm (9 1/2")





Ligature Needles
Unterbindungsnadeln
Agujas para ligaduras
Aiguilles à ligatures
Aghi per legatura

Deschamps Ligature Needle
Blunt For Left Hand

Catalog # Size
SU-873-24 24 cm (9 1/2")



Deschamps Ligature Needle
Blunt For Right Hand

Catalog # Size
SU-875-20 20 cm (8")



Deschamps Ligature Needle
Blunt For Left Hand

Catalog # Size
SU-877-20 20 cm (8")



Deschamps Ligature Needle
Sharp For Right Hand

Catalog # Size
SU-879-20 20 cm (8")



Deschamps Ligature Needle
Sharp For Left Hand

Catalog # Size
SU-881-20 20 cm (8")





Ligature Needles
Unterbindungsnadeln
Agujas para ligaduras
Aiguilles à ligatures
Aghi per legatura

Deschamps Ligature Needle
Blunt For Right Hand

Catalog #	Size
SU-883-21	20 cm (8")



Deschamps Ligature Needle
Blunt For Left Hand

Catalog #	Size
SU-885-21	20 cm (8")



Deschamps Ligature Needle
Sharp For Right Hand

Catalog #	Size
SU-887-20	20 cm (8")



Deschamps Ligature Needle
Sharp For Left Hand

Catalog #	Size
SU-889-20	20 cm (8")



Deschamps Ligature Needle
Sharp For Right Hand

Catalog #	Size
SU-891-21	20 cm (8")



INCHES
MILLIMETERS



Ligature Needles
Unterbindungsnadeln
Agujas para ligaduras
Aiguilles à ligatures
Aghi per legatura

Deschamps Ligature Needle

Sharp For Left Hand



Catalog # Size
SU-893-21 20 cm (8")

Deschamps Ligature Needle

Blunt For Right Hand



Catalog # Size
SU-895-27 27 cm (10 3/4")

Deschamps Ligature Needle

Blunt For Left Hand



Catalog # Size
SU-897-27 27 cm (10 3/4")

Deschamps Ligature Needle

Sharp For Right Hand



Catalog # Size
SU-901-27 27 cm (10 3/4")

Deschamps Ligature Needle

Sharp For Left Hand



Catalog # Size
SU-903-27 27 cm (10 3/4")



Ligature Needles
Unterbindungsnadeln
Agujas para ligaduras
Aiguilles à ligatures
Aghi per legatura

Deschamps Ligature Needle
Blunt For Right Hand



Catalog #	Size
SU-905-28	28 cm (11")

Deschamps Ligature Needle
Blunt For Left Hand



Catalog #	Size
SU-907-28	28 cm (11")

Deschamps Ligature Needle
Sharp For Right Hand



Catalog #	Size
SU-911-20	20 cm (8")

Deschamps Ligature Needle
Sharp For Left Hand



Catalog #	Size
SU-913-20	20 cm (8")

Deschamps Ligature Needle
Blunt For Right Hand




Catalog #	Size
SU-915-20	20 cm (8")

INCHES
MILLIMETERS



Ligature Needles
Unterbindungsnadeln
Agujas para ligaduras
Aiguilles à ligatures
Aghi per legatura

Deschamps Ligature Needle
Blunt For Left Hand



Catalog # Size
SU-917-20 20 cm (8")

Deschamps Ligature Needle
Blunt For Right Hand




Catalog # Size
SU-921-20 21 cm (8 1/4")

Deschamps Ligature Needle
Blunt For Left Hand




Catalog # Size
SU-923-20 21 cm (8 1/4")

Hurd Ligature Needle
Sharp For Right Hand



Catalog # Size
SU-927-21 20.5 cm (8")

Hurd Ligature Needle
Sharp For Left Hand




Catalog # Size
SU-929-21 20.5 cm (8")



Ligature Needles
Unterbindungsnadeln
Agujas para ligaduras
Aiguilles à ligatures
Aghi per legatura

Brunner Ligature Needle



Catalog #	Size
SU-931-01	30 cm (11 3/4")
SU-931-02	30 cm (11 3/4")
SU-931-03	30 cm (11 3/4")
SU-931-04	30 cm (11 3/4")




Kirschner Ligature Needle



Catalog #	Size
SU-935-25	25.5 cm (10")




Ligature Needle



Catalog #	Size
SU-937-28	28 cm (11")



Martin Ligature Needle



Catalog #	Size
SU-939-25	19.5 cm (7 3/4")



INCHES
MILLIMETERS



Ligature Conductors
Führungshohlsonden
Sondas conductoras
Sondes conductrices
Sonde conduttrici

Brunner Ligature Conductor
Ferrozell



Catalog #	Size
SU-941-18	18 cm (7")

Kirschner Ligature Conductor
Ferrozell



Catalog #	Size
SU-943-27	27 cm (10 3/4")

Ligature Conductor
Ferrozell




Catalog #	Size
SU-945-19	19 cm (7 1/2")

Ligature Conductor
Ferrozell



Catalog #	Size
SU-949-24	24 cm (9 1/2")

Schmieden Ligature Conductor




Catalog #	Size
SU-951-19	18 cm (7")




Ligature Conductors
Führungshohlsonden
Sondas conductoras
Sondes conductrices
Sonde conduttrici

Schmieden-Payr Ligature Conductor



Catalog #	Size
SU-953-19	18 cm (7")


Koeing Ligature Conductor



Catalog #	Size
SU-955-03	19.5 cm (7 3/4")
SU-955-05	19.5 cm (7 3/4")
SU-955-07	19.5 cm (7 3/4")




Payr Ligature Conductor



Catalog #	Size
SU-957-22	22 cm (8 3/4")



Payr Ligature Conductor



Catalog #	Size
SU-959-22	22 cm (8 3/4")





Modell Wien Ligature Catcher



Catalog #	Size
SU-961-23	23.5 cm (9 1/4")

Bozemann Knot Tier



Catalog #	Size
SU-963-21	21 cm (8 1/4")

Ligature Adjuster



Catalog #	Size
SU-965-23	23 cm (9")

Caspar Ligature Adjuster
Bayonet Shaped



Catalog #	Size
SU-967-23	23 cm (9")

Millin Ligature Holding Forcep



Catalog #	Size
SU-969-24	24 cm (9 1/2")



Clip Applying Forceps
Wundklammerpinzetten, Wundklammerzängchen
Pinzas aplicadoras para agrafes
Pinces à poser les agrafes
Pinze per applicare le agrafe

Childe Approximation Forcep
Heavy Pattern With Clip Holder



Catalog #	Size
SU-980-18	17.5 cm (7")

Michel Clip Applying Forcep



Catalog #	Size
SU-984-12	12 cm (4 3/4")

Hegenbarth Clip Applying Forcep



Catalog #	Size
SU-986-12	12 cm (4 3/4")

Wachenfeldt Clip Applying Forcep



Catalog #	Size
SU-988-12	12 cm (4 3/4")

Michel Clip Applying & Removing Forcep
For Applying and Removal of Clips



Catalog #	Size
SU-990-12	12.5 cm (5")



Clip Applying Forceps
Wundklammerpinzetten, Wundklammerzängchen
Pinzas aplicadoras para agrafes
Pincas à poser les agrafes
Pinze per applicare le agrafe

Michel Clip Applying & Removing Forcep
Long Jaws For Applying and Removal of Clips



Catalog #	Size
SU-992-12	13 cm (5")

Collin Clip Applying Forcep
With Screw Lock



Catalog #	Size
SU-994-13	13 cm (5")

Richter Clip Applying Forcep
With Screw Lock



Catalog #	Size
SU-997-14	13.5 cm (5 1/4")

Michel Suture Clips
Pack of 1000

Catalog #	Clip Size
SU-999-01	7.5 x 1.75 mm
SU-999-02	11 x 2.0 mm
SU-999-03	12 x 3.0 mm
SU-999-04	14 x 3.0 mm
SU-999-05	16 x 3.0 mm
SU-999-06	18 x 3.0 mm
SU-999-07	20 x 3.0 mm
SU-999-08	22 x 3.0 mm

