



SURTEX
INSTRUMENTS

Crafted to Excellence



SS

SCISSORS

Surtex scissors are renowned for their **Xedge**™ technology enabling users to experience a new dimension in cutting and dissecting performance. Surtex unique **Xedge**™ leading cutting technology allows scissors to self sharpen as the design ensures the scissor hinge provides spring tension on the two blades so, as they cut material, one blade presses against the other, sharpening the surface as they close. Surtex **Xedge**™ scissors permit minimum force but confer maximum cutting strength in all types of soft tissues even after pro-longed use.

In addition to standard **Xedge**™ our scissors are available in the following finishes:

SUPER Edge™ Super Cut Scissors

These scissors deliver unbeatable performance when cutting soft tissues because one blade is manufactured with micro-serrations perpendicular to the cutting edge and the other blade has a bevel-honed edge resulting in scissors that cut forward.

The special micro serrations in conjunction with the uniquely ground edges ensure precise, fine and minimally traumatizing cuts. All popular models are available in this finish and identified by black ring handles.

ULTRA Cut™ Tungsten Carbide Scissors

These scissors feature carbide inserts a material tougher than stainless steel enabling the scissors to last up to 6 times longer than standard inserts. These scissors are recommended for use in the cutting of cartilage and fur. All popular models are available in this finish and identified by 18K gold plated ring handles.

XTS Cut™ Super Cut + Tungsten Carbide Scissors

These are the best scissors in our range and combine features of both **ULTRA Cut**™ & **SUPER Edge**™ featuring a combination of tungsten carbide inserts and special micro serrations to deliver the ultimate in cutting performance. All popular model are available in this finish and identified by 18K gold plated & black ring handles.

WARNING: **XTS Cut**™ Super Cut Scissors must not be used for cutting wire or other hard objects. They should also not be used for cutting hair, fur, or tendons as this may irreversibly damage the blades and invalidate your warranty.

Scheren

Tijeras

Torbici

Ciseaux



Operating Scissors
Chirurgische Scheren
Tijeras quirúrgicas
Ciseaux chirurgicaux
Forbici per chirurgia


Operating Scissors
Straight, Blunt/Blunt

Catalog #	Size
SS-002-11	12 cm (4 3/4")
SS-002-13	13 cm (5")
SS-002-14	14.5 cm (5 3/4")
SS-002-16	16.5 cm (6 1/2")
SS-002-18	18.5 cm (7 1/4")
SS-002-20	20.5 cm (8")




Operating Scissors
Curved, Blunt/Blunt

Catalog #	Size
SS-003-11	12 cm (4 3/4")
SS-003-13	13 cm (5")
SS-003-14	14.5 cm (5 3/4")
SS-003-16	16.5 cm (6 1/2")
SS-003-18	18.5 cm (7 1/4")
SS-003-20	20.5 cm (8")



Operating Scissors
Straight, Sharp/Blunt

Catalog #	Size
SS-004-11	12 cm (4 3/4")
SS-004-13	13 cm (5")
SS-004-14	14.5 cm (5 3/4")
SS-004-16	16.5 cm (6 1/2")
SS-004-18	18.5 cm (7 1/4")
SS-004-20	20.5 cm (8")





Operating Scissors
Chirurgische Scheren
Tijeras quirúrgicas
Ciseaux chirurgicaux
Forbici per chirurgia

SS

SCISSORS

Operating Scissors
Straight, Sharp/Blunt


SUPER*Edge™*

Catalog #	Size
SS-004-13 SC	13 cm (5")
SS-004-14 SC	14.5 cm (5 ¾")



Operating Scissors
Curved, Sharp/Blunt

Catalog #	Size
SS-005-11	12 cm (4 ¾")
SS-005-13	13 cm (5")
SS-005-14	14.5 cm (5 ¾")
SS-005-16	16.5 cm (6 ½")
SS-005-18	18.5 cm (7 ¼")
SS-005-20	20.5 cm (8")



Operating Scissors
Curved, Sharp/Blunt,

SUPER*Edge™*

Catalog #	Size
SS-005-13 SC	13 cm (5")
SS-005-14 SC	14.5 cm (5 ¾")



Operating Scissors
Straight, Sharp/Sharp

Catalog #	Size
SS-006-11	12 cm (4 ¾")
SS-006-13	13 cm (5")
SS-006-14	14.5 cm (5 ¾")
SS-006-16	16.5 cm (6 ½")
SS-006-18	18.5 cm (7 ¼")
SS-006-20	20.5 cm (8")






Operating Scissors
Chirurgische Scheren
Tijeras quirúrgicas
Ciseaux chirurgicaux
Forbici per chirurgia

Operating Scissors
Curved, Sharp/Sharp

Catalog #	Size
SS-007-11	12 cm (4 ¾")
SS-007-13	13 cm (5")
SS-007-14	14.5 cm (5 ¾")
SS-007-16	16.5 cm (6 ½")
SS-007-18	18.5 cm (7 ¼")
SS-007-20	20.5 cm (8")




Operating Scissors
Straight, Blunt/Blunt

Catalog #	Size
SS-010-14	14.5 cm (5 ¾")




Operating Scissors
Curved, Blunt/Blunt

Catalog #	Size
SS-011-14	14.5 cm (5 ¾")



Operating Scissors
Straight, Sharp/Blunt

Catalog #	Size
SS-012-14	14.5 cm (5 ¾")






Operating Scissors
Chirurgische Scheren
Tijeras quirúrgicas
Ciseaux chirurgicaux
Forbici per chirurgia

SS

SCISSORS


Operating Scissors
Curved, Sharp/Blunt

Catalog #	Size
SS-013-14	14.5 cm (5 ¾")




Operating Scissors
Straight, Sharp/Sharp

Catalog #	Size
SS-014-14	14.5 cm (5 ¾")



Operating Scissors
Curved, Sharp/Sharp

Catalog #	Size
SS-015-14	14.5 cm (5 ¾")



Incision Scissor
Straight

Catalog #	Size
SS-018-13	13 cm (5")
SS-020-14	14.5 cm (5 ¾")
SS-022-16	16.5 cm (6 ½")



Systrunk Nail Splitting Scissor
Straight


Catalog #	Size
SS-024-13	13 cm (5")





Operating Scissors, Dissecting Scissors
Chirurgische Scheren, Präparierscheren
Tijeras quirúrgicas, Tijeras de disección
Ciseaux chirurgicaux, Ciseaux à dissection
Forbici per chirurgia, Forbici per dissezione

Deaver Operating Scissor
Straight, Sharp/Blunt




Catalog #	Size
SS-030-14	14 cm (5 1/2")

Deaver Operating Scissor
Curved, Sharp/Blunt



Catalog #	Size
SS-031-14	14 cm (5 1/2")

Mayo Dissecting Scissor
Straight



Catalog #	Size
SS-040-14	14.5 cm (5 3/4")
SS-040-15	15 cm (6")
SS-040-17	17 cm (6 3/4")
SS-040-23	23.5 cm (9 1/4")

Mayo Dissecting Scissor
Straight



SUPER Edge™

Catalog #	Size
SS-040-14 SC	14.5 cm (5 3/4")
SS-040-15 SC	15 cm (6")
SS-040-17 SC	17 cm (6 3/4")




Dissecting Scissors
Präparierscheren
Tijeras de disección
Ciseaux à dissection
Forbici per dissezione

SS

SCISSORS

Mayo Dissecting Scissor
Curved


Catalog #	Size
SS-041-14	14.5 cm (5 ¾")
SS-041-15	15 cm (6")
SS-041-17	17 cm (6 ¾")
SS-041-23	23.5 cm (9 ¼")



Mayo Dissecting Scissor
Curved

SUPER *Edge™*

Catalog #	Size
SS-041-14 SC	14.5 cm (5 ¾")
SS-041-15 SC	15 cm (6")
SS-041-17 SC	17 cm (6 ¾")




Mayo Dissecting Scissor
Straight, With Chamfered Blades

Catalog #	Size
SS-042-19	19 cm (7 ½")
SS-042-21	21.5 cm (8 ½")



Mayo Dissecting Scissor
Curved, With Chamfered Blades

Catalog #	Size
SS-043-19	19 cm (7 ½")
SS-043-21	21.5 cm (8 ½")





Dissecting Scissors
Präparierscheren
Tijeras de disección
Ciseaux à dissection
Forbici per dissezione

Mayo-Harrington Dissecting Scissor
Straight



Catalog #	Size
SS-046-23	23 cm (9")

Mayo-Harrington Dissecting Scissor
Straight



SUPER Edge™


Catalog #	Size
SS-046-23 SC	23 cm (9")

Mayo-Harrington Dissecting Scissor
Curved



Catalog #	Size
SS-047-23	22.5 cm (9")
SS-047-28	28 cm (11")

Mayo-Harrington Dissecting Scissor
Curved



SUPER Edge™

Catalog #	Size
SS-047-23 SC	22.5 cm (9")



Dissecting Scissors
Präparierscheren
Tijeras de disección
Ciseaux à dissection
Forbici per dissezione

SS

SCISSORS

Mayo-Stille Dissecting Scissor
Straight

Catalog #	Size
SS-050-15	15 cm (6")
SS-050-17	17 cm (6 ¾")
SS-050-19	19.5 cm (7 ¾")
SS-050-21	21.5 cm (8 ½")



Mayo-Stille Dissecting Scissor
Straight


SUPER Edge™

Catalog #	Size
SS-050-15 SC	15 cm (6")
SS-050-17 SC	17 cm (6 ¾")



Mayo-Stille Dissecting Scissor
Curved

Catalog #	Size
SS-051-15	15 cm (6")
SS-051-17	17 cm (6 ¾")
SS-051-19	19.5 cm (7 ¾")
SS-051-21	21.5 cm (8 ½")



Mayo-Stille Dissecting Scissor
Curved

SUPER Edge™

Catalog #	Size
SS-051-15 SC	15 cm (6")
SS-051-17 SC	17 cm (6 ¾")





Dissecting Scissors, Gynecological Scissors
Präparierscheren, Gynäkologische Scheren
Tijeras de disección, Tijeras para ginecología
Ciseaux à dissection, Ciseaux à gynécologie
Forbici per dissezione, Forbici per ginecologia

Mayo-Stille Dissecting Scissor
Straight, With Chamfered Blades

Catalog #	Size
SS-052-14	14 cm (5 1/2")
SS-052-17	17 cm (6 3/4")



Mayo-Stille Dissecting Scissor
Curved, With Chamfered Blades

Catalog #	Size
SS-053-14	14 cm (5 1/2")
SS-053-17	17 cm (6 3/4")




Stille Dissecting Scissor
Straight

Catalog #	Size
SS-056-15	15 cm (6")



Stille Dissecting Scissor
Curved

Catalog #	Size
SS-057-15	15 cm (6")



Mayo-Noble Gynecological Scissor
Straight

Catalog #	Size
SS-060-17	16.5 cm (6 1/2")





Gynecological Scissors
Gynäkologische Scheren
Tijeras para ginecología
Ciseaux à gynécologie
Forbici per ginecologia

SS

SCISSORS

Mayo-Noble Gynecological Scissor
Curved

Catalog #	Size
SS-061-17	16.5 cm (6 1/2")



Sims Gynecological Scissor
Straight

Catalog #	Size
SS-070-20	20 cm (8")
SS-070-23	23 cm (9")



Sims Gynecological Scissor
Straight

SUPER Edge™

Catalog #	Size
SS-070-20 SC	20 cm (8")
SS-070-23 SC	23 cm (9")



Sims Gynecological Scissor
Curved

Catalog #	Size
SS-071-20	20 cm (8")
SS-071-23	23 cm (9")





Gynecological Scissors
Gynäkologische Scheren
Tijeras para ginecología
Ciseaux à gynécologie
Forbici per ginecologia

Sims Gynecological Scissor
Curved


SUPER Edge™

Catalog #	Size
SS-071-20 SC	20 cm (8")
SS-071-23 SC	23 cm (9")



Doyen Gynecological Scissor
Straight

Catalog #	Size
SS-074-17	18.5 cm (7 1/4")



Doyen Gynecological Scissor
Curved

Catalog #	Size
SS-075-17	18.5 cm (7 1/4")



Sims-Siebold Gynecological Scissor
Straight


Catalog #	Size
SS-078-24	24.5 cm (9 3/4")



Sims-Siebold Gynecological Scissor
Straight

SUPER Edge™

Catalog #	Size
SS-078-24 SC	24.5 cm (9 3/4")






Gynecological Scissors
Gynäkologische Scheren
Tijeras para ginecología
Ciseaux à gynécologie
Forbici per ginecologia

SS

SCISSORS

Sims-Siebold Gynecological Scissor
Curved


Catalog #	Size
SS-079-24	24.5 cm (9 3/4")



Sims-Siebold Gynecological Scissor
Curved

SUPER *Edge*[™]

Catalog #	Size
SS-079-24 SC	24.5 cm (9 3/4")




Sims Gynecological Scissor
Straight, Sharp/Sharp

Catalog #	Size
SS-082-20	20 cm (8")




Sims Gynecological Scissor
Curved, Blunt/Blunt

Catalog #	Size
SS-085-21	21 cm (8 1/4")



Wertheim Gynecological Scissor
Curved

Catalog #	Size
SS-087-20	20 cm (8")
SS-087-23	23 cm (9")






Gynecological Scissors, Dissecting Scissors
Gynäkologische Scheren, Präparierscheren
Tijeras para ginecología, Tijeras de disección
Ciseaux à gynécologie, Ciseaux à dissection
Forbici per ginecologia, Forbici per dissezione


Bozemann Gynecological Scissor
S Shaped

Catalog #	Size
SS-089-22	22 cm (8 ¾")




Kieback Parametrium Scissor
Angled 90°

Catalog #	Size
SS-091-24	24 cm (9 ½")
SS-091-35	35 cm (13 ¾")




Kieback Parametrium Scissor
Curved to Side, One Probe Pointed Blade
One Toothed Cutting Edge

Catalog #	Size
SS-093-26	26 cm (10 ¼")
SS-093-35	35 cm (13 ¾")



Kelly Dissecting Scissor
Straight


Catalog #	Size
SS-100-16	16 cm (6 ¼")



Kelly Dissecting Scissor
Straight

SUPER Edge™

Catalog #	Size
SS-100-16 SC	16 cm (6 ¼")






Dissecting Scissors
Präparierscheren
Tijeras de disección
Ciseaux à dissection
Forbici per dissezione

SS

SCISSORS

Kelly Dissecting Scissor
Curved


Catalog #	Size
SS-101-16	16 cm (6 1/4")



Kelly Dissecting Scissor
Curved


SUPER *Edge*[™]

Catalog #	Size
SS-101-16 SC	16 cm (6 1/4")



Joseph Dissecting Scissor
Straight


Catalog #	Size
SS-104-15	15 cm (6")



Joseph Dissecting Scissor
Straight


SUPER *Edge*[™]

Catalog #	Size
SS-104-15 SC	15 cm (6")



Joseph Dissecting Scissor
Curved

Catalog #	Size
SS-105-15	15 cm (6")






Dissecting Scissors
Präparierscheren
Tijeras de disección
Ciseaux à dissection
Forbici per dissezione

Joseph Dissecting Scissor
Curved

SUPER Edge™

Catalog #	Size
SS-105-15 SC	15 cm (6")




Weller Dissecting Scissor
Heavy Pattern

Catalog #	Size
SS-111-28	28 cm (11")




Resano Cartilage Scissor

Catalog #	Size
SS-113-25	25 cm (9 3/4")



Metzenbaum Dissecting Scissor
Straight, Blunt/Blunt


Catalog #	Size
SS-120-11	11.5 cm (4 1/2")



Metzenbaum Dissecting Scissor
Straight, Blunt/Blunt

SUPER Edge™

Catalog #	Size
SS-120-11 SC	11.5 cm (4 1/2")





Dissecting Scissors
Präparierscheren
Tijeras de disección
Ciseaux à dissection
Forbici per dissezione

SS

SCISSORS

Metzenbaum Dissecting Scissor
Straight, Blunt/Blunt


SUPER*Edge*[™]

Catalog #	Size
SS-122-15 SC	15.5 cm (6")



Metzenbaum Dissecting Scissor
Curved, Blunt/Blunt

Catalog #	Size
SS-123-15	15.5 cm (6")



Metzenbaum Dissecting Scissor
Curved, Blunt/Blunt

SUPER*Edge*[™]

Catalog #	Size
SS-123-15 SC	15.5 cm (6")



Metzenbaum Dissecting Scissor
Straight, Blunt/Blunt

Catalog #	Size
SS-130-14	14.5 cm (5 ¾")



Metzenbaum Dissecting Scissor
Straight, Blunt/Blunt

SUPER*Edge*[™]

Catalog #	Size
SS-130-14 SC	14.5 cm (5 ¾")






Dissecting Scissors
Präparierscheren
Tijeras de disección
Ciseaux à dissection
Forbici per dissezione

Metzenbaum Dissecting Scissor
Curved, Blunt/Blunt


Catalog #	Size
SS-121-11	11.5 cm (4 1/2")



Metzenbaum Dissecting Scissor
Curved, Blunt/Blunt

SUPER Edge™

Catalog #	Size
SS-121-11 SC	11.5 cm (4 1/2")



Metzenbaum Dissecting Scissor
Straight, Blunt/Blunt

Catalog #	Size
SS-122-15	15.5 cm (6")



Metzenbaum Dissecting Scissor
Curved, Blunt/Blunt

Catalog #	Size
SS-131-14	14.5 cm (5 3/4")



Metzenbaum Dissecting Scissor
Curved, Blunt/Blunt

SUPER Edge™

Catalog #	Size
SS-131-14 SC	14.5 cm (5 3/4")






Dissecting Scissors
Präparierscheren
Tijeras de disección
Ciseaux à dissection
Forbici per dissezione

SS

SCISSORS


Metzenbaum Dissecting Scissor
Straight, Sharp/Blunt

Catalog #	Size
SS-132-14	14.5 cm (5 ¾")



Metzenbaum Dissecting Scissor
Curved, Sharp/Blunt

Catalog #	Size
SS-133-14	14.5 cm (5 ¾")



Metzenbaum Dissecting Scissor
Straight, Sharp/Sharp

Catalog #	Size
SS-134-14	14.5 cm (5 ¾")



Metzenbaum Dissecting Scissor
Curved, Sharp/Sharp

Catalog #	Size
SS-135-14	14.5 cm (5 ¾")






Dissecting Scissors
Präparierscheren
Tijeras de disección
Ciseaux à dissection
Forbici per dissezione

Metzenbaum-Nelson Dissecting Scissor
Straight, Blunt/Blunt


Catalog #	Size
SS-140-18	18 cm (7")
SS-140-20	20.5 cm (8")
SS-140-23	23 cm (9")
SS-140-26	26 cm (10 1/4")
SS-140-28	28.5 cm (11 1/4")
SS-140-30	31 cm (12 1/4")



Metzenbaum-Nelson Dissecting Scissor
Straight, Blunt/Blunt



SUPER Edge™

Catalog #	Size
SS-140-18 SC	18 cm (7")
SS-140-20 SC	20.5 cm (8")
SS-140-23 SC	23 cm (9")
SS-140-26 SC	26 cm (10 1/4")



Metzenbaum-Nelson Dissecting Scissor
Curved, Blunt/Blunt

Catalog #	Size
SS-141-18	18 cm (7")
SS-141-20	20.5 cm (8")
SS-141-23	23 cm (9")
SS-141-26	26 cm (10 1/4")
SS-141-28	28.5 cm (11 1/4")
SS-141-30	31 cm (12 1/4")





Dissecting Scissors
Präparierscheren
Tijeras de disección
Ciseaux à dissection
Forbici per dissezione


SS

SCISSORS

Metzenbaum-Nelson Dissecting Scissor
Curved, Blunt/Blunt


SUPER *Edge*[™]

Catalog #	Size
SS-141-18 SC	18 cm (7")
SS-141-20 SC	20.5 cm (8")
SS-141-23 SC	23 cm (9")
SS-141-26 SC	26 cm (10 ¼")




Metzenbaum-Nelson Dissecting Scissor
Straight, Sharp/Sharp

Catalog #	Size
SS-144-18	18 cm (7")
SS-144-20	20.5 cm (8")




Metzenbaum Dissecting Scissor
Straight, Sharp/Sharp, Toothed

Catalog #	Size
SS-146-14	14.5 cm (5 ¾")



Metzenbaum Dissecting Scissor
Curved, Sharp/Sharp, Toothed

Catalog #	Size
SS-147-14	14.5 cm (5 ¾")






Dissecting Scissors, Face-lift Scissors
Präparierscheren, Face-lift-Scheren
Tijeras de disección, Tijeras para Face-Lift
Ciseaux à dissection, Ciseaux pour Face-lift
Forbici per dissezione, Forbici per Face-lift


Metzenbaum Nerve Dissecting Scissor
Straight

Catalog #	Size
SS-148-15	15.5 cm (6")




Metzenbaum Nerve Dissecting Scissor
Curved

Catalog #	Size
SS-149-15	15.5 cm (6")




McIndoe Cartilage Scissor
Straight, Toothed

Catalog #	Size
SS-154-18	18.5 cm (7 1/4")



McIndoe Cartilage Scissor
Curved, Toothed

Catalog #	Size
SS-155-18	18.5 cm (7 1/4")



Kaye Face-Lift Scissor
Curved

SUPER *Edge*[™]

Catalog #	Size
SS-161-18 SC	18.5 cm (7 1/4")





Face-lift Scissors, Dissecting Scissors
Face-lift-Scheren, Präparierscheren
Tijeras para Face-Lift, Tijeras de disección
Ciseaux pour Face-lift, Ciseaux à dissection
Forbici per Face-lift, Forbici per dissezione


SS

SCISSORS

Gorney Face-Lift Scissor
Straight

SUPER Edge™

Catalog #	Size
SS-166-12 SC	12.5 cm (5")
SS-166-19 SC	19.5 cm (7 ¾")



Gorney Face-Lift Scissor
Curved

SUPER Edge™

Catalog #	Size
SS-167-19 SC	19.5 cm (7 ¾")



Aston Face-Lift Scissor
Straight

SUPER Edge™

Catalog #	Size
SS-170-17 SC	17 cm (6 ¾")



Lexer Dissecting Scissor
Straight

Catalog #	Size
SS-176-16	16 cm (6 ¼")



Lexer Dissecting Scissor
Straight

SUPER Edge™

Catalog #	Size
SS-176-16 SC	16 cm (6 ¼")





Dissecting Scissors, Nasal Scissors
Präparierscheren, Nasenscheren
Tijeras de disección, Tijeras nasales
Ciseaux à dissection, Ciseaux nasaux
Forbici per dissezione, Forbici nasali

Lexer Dissecting Scissor
Curved

Catalog #	Size
SS-177-16	16 cm (6 1/4")



Lexer Dissecting Scissor
Curved


SUPER Edge™

Catalog #	Size
SS-177-16 SC	16 cm (6 1/4")




Cottle-Masing Plastic Surgery Scissor
Curved, Sharp/Sharp

Catalog #	Size
SS-181-10	10.5 cm (4 1/4")



Cottle-Masing Plastic Surgery Scissor
Curved, Blunt/Blunt

Catalog #	Size
SS-183-10	10.5 cm (4 1/4")



Cinelli Nasal Scissor
Curved

Catalog #	Size
SS-187-11	11.5 cm (4 1/2")





Nasal Scissors, Delicate Dissecting Scissors
Nasenscheren, Feine Präparierscheren
Tijeras nasales, Tijeras finas de disección
Ciseaux nasaux, Ciseaux fins à dissection
Forbici nasali, Forbici fini per dissezione

Fomon Nasal Scissor
Curved

Catalog #	Size
SS-189-14	14 cm (5 ½")




Thorek Dissecting Scissor
Curved

Catalog #	Size
SS-191-19	18.5 cm (7 ¼")




Lexer-Baby Dissecting Scissor
Straight

Catalog #	Size
SS-194-10	10 cm (4")




Lexer-Baby Dissecting Scissor
Curved

Catalog #	Size
SS-195-10	10 cm (4")



Lexer-Fine Delicate Dissecting Scissor
Straight, Slender Pattern

Catalog #	Size
SS-196-16	16.5 cm (6 ½")





Delicate Dissecting Scissors
Feine Präparierscheren
Tijeras finas de disección
Ciseaux fins à dissection
Forbici fini per dissezione

Lexer-Fine Delicate Dissecting Scissor
Curved, Slender Pattern



Catalog #	Size
SS-197-16	16.5 cm (6 1/2")

Metzenbaum-Fine Delicate Dissecting Scissor
Straight, Blunt/Blunt, Slender Pattern



Catalog #	Size
SS-200-14	14.5 cm (5 3/4")


Metzenbaum-Fine Delicate Dissecting Scissor
Straight, Blunt/Blunt, Slender Pattern



SUPER Edge™

Catalog #	Size
SS-200-14 SC	14.5 cm (5 3/4")

Metzenbaum-Fine Delicate Dissecting Scissor
Curved, Blunt/Blunt, Slender Pattern



Catalog #	Size
SS-201-14	14.5 cm (5 3/4")

Metzenbaum-Fine Delicate Dissecting Scissor
Curved, Blunt/Blunt, Slender Pattern



SUPER Edge™

Catalog #	Size
SS-201-14 SC	14.5 cm (5 3/4")




Delicate Dissecting Scissors
Feine Präparierscheren
Tijeras finas de disección
Ciseaux fins à dissection
Forbici fini per dissezione

SS

SCISSORS


Metzenbaum-Fine Delicate Dissecting Scissor
Straight, Sharp/Sharp, Slender Pattern

Catalog #	Size
SS-202-14	14.5 cm (5 ¾")




Metzenbaum-Fine Delicate Dissecting Scissor
Curved, Sharp/Sharp, Slender Pattern

Catalog #	Size
SS-203-14	14.5 cm (5 ¾")



Metzenbaum-Fine Delicate Dissecting Scissor
Straight, Blunt/Blunt, Slender Pattern


Catalog #	Size
SS-210-18	18 cm (7")
SS-210-20	20 cm (8")
SS-210-23	23 cm (9")



Metzenbaum-Fine Delicate Dissecting Scissor
Straight, Blunt/Blunt, Slender Pattern

SUPER *Edge*™

Catalog #	Size
SS-210-18 SC	18 cm (7")
SS-210-20 SC	20 cm (8")
SS-210-23 SC	23 cm (9")






Delicate Dissecting Scissors, Dissecting Scissors
Feine Präparierscheren, Präparierscheren
Tijeras finas de disección, Tijeras de disección
Ciseaux fins à dissection, Ciseaux à dissection
Forbici fini per dissezione, Forbici per dissezione

Metzenbaum-Fine Delicate Dissecting Scissor
Curved, Blunt/Blunt, Slender Pattern


Catalog #	Size
SS-211-18	18 cm (7")
SS-211-20	20 cm (8")
SS-211-23	23 cm (9")



Metzenbaum-Fine Delicate Dissecting Scissor
Curved, Blunt/Blunt, Slender Pattern


SUPER Edge™

Catalog #	Size
SS-211-18 SC	18 cm (7")
SS-211-20 SC	20 cm (8")
SS-211-23 SC	23 cm (9")



Metzenbaum-Fine Delicate Dissecting Scissor
Straight, Sharp/Sharp, Slender Pattern

Catalog #	Size
SS-214-18	18 cm (7")
SS-214-20	20 cm (8")
SS-214-23	23 cm (9")



Metzenbaum-Fine Delicate Dissecting Scissor
Straight, Sharp/Sharp, Slender Pattern

SUPER Edge™

Catalog #	Size
SS-214-18 SC	18 cm (7")
SS-214-20 SC	20 cm (8")
SS-214-23 SC	23 cm (9")





Delicate Dissecting Scissors, Dissecting Scissors
Feine Präparierscheren, Präparierscheren
Tijeras finas de disección, Tijeras de disección
Ciseaux fins à dissection, Ciseaux à dissection
Forbici fini per dissezione, Forbici per dissezione

SS

SCISSORS

Metzenbaum-Fine Delicate Dissecting Scissor
Curved, Sharp/Sharp, Slender Pattern


Catalog #	Size
SS-215-18	18 cm (7")
SS-215-20	20 cm (8")
SS-215-23	23 cm (9")



Metzenbaum-Fine Delicate Dissecting Scissor
Curved, Sharp/Sharp, Slender Pattern


SUPER Edge™

Catalog #	Size
SS-215-18 SC	18 cm (7")
SS-215-20 SC	20 cm (8")
SS-215-23 SC	23 cm (9")




Metzenbaum-Fine Delicate Dissecting Scissor
Curved, S Shaped

Catalog #	Size
SS-219-18	18 cm (7")
SS-219-20	20 cm (8")
SS-219-23	23 cm (9")



Salzer Cleft Palate Dissecting Scissor
Curved, Fine Pattern

Catalog #	Size
SS-223-12	12.5 cm (5")
SS-223-18	18 cm (7")






Dissecting Scissors, Face-lift Scissors
Präparierscheren, Face-lift-Scheren
Tijeras de disección, Tijeras para Face-Lift
Ciseaux à dissection, Ciseaux pour Face-lift
Forbici per dissezione, Forbici per Face-lift


Sanvenero Dissecting Scissor
Straight

Catalog #	Size
SS-226-14	14 cm (5 ½")




Sanvenero Dissecting Scissor
Curved

Catalog #	Size
SS-227-14	14 cm (5 ½")




Tenotomy Scissor
Straight

Catalog #	Size
SS-230-11	11.5 cm (4 ½")




Tenotomy Scissor
Curved

Catalog #	Size
SS-231-11	11.5 cm (4 ½")



Peck-Joseph Face-Lift Scissor
Straight, Sharp/Sharp

Catalog #	Size
SS-236-14	14.5 cm (5 ¾")







Face-lift Scissors, Dissecting Scissors
Face-lift-Scheren, Präparierscheren
Tijeras para Face-Lift, Tijeras de disección
Ciseaux pour Face-lift, Ciseaux à dissection
Forbici per Face-lift, Forbici per dissezione

SS



SCISSORS

Peck-Joseph Face-Lift Scissor
Curved, Sharp/Sharp




Catalog #	Size
SS-237-14	14.5 cm (5 ¾")

Peck-Joseph Face-Lift Scissor
Straight, Blunt/Blunt





Catalog #	Size
SS-238-14	14.5 cm (5 ¾")

Peck-Joseph Face-Lift Scissor
Curved, Blunt/Blunt




Catalog #	Size
SS-239-14	14.5 cm (5 ¾")

Shea Vein Graft Scissor
Curved



Catalog #	Size
SS-241-12	12 cm (4 ¾")

Ragnell (Kilner) Dissecting Scissor
Curved




Catalog #	Size
SS-243-12	12 cm (4 ¾")
SS-243-15	15 cm (6")
SS-243-18	18 cm (7")



Dissecting Scissors
Präparierscheren
Tijeras de disección
Ciseaux à dissection
Forbici per dissezione

Nerve Dissecting Scissor
Curved, Semi Sharp

Catalog #	Size
SS-245-15	15 cm (6")



Dissecting Scissor
Straight

Catalog #	Size
SS-248-14	14.5 cm (5 3/4")




Dissecting Scissor
Curved

Catalog #	Size
SS-249-14	14.5 cm (5 3/4")



Reynolds Dissecting Scissor
Straight


Catalog #	Size
SS-250-15	15.5 cm (6")
SS-250-18	18 cm (7")



Reynolds Dissecting Scissor
Straight

SUPER Edge™

Catalog #	Size
SS-250-15 SC	15.5 cm (6")
SS-250-18 SC	18 cm (7")






Dissecting Scissors
Präparierscheren
Tijeras de disección
Ciseaux à dissection
Forbici per dissezione

SS


SCISSORS

Reynolds Dissecting Scissor
Curved




Catalog #	Size
SS-251-15	15.5 cm (6")
SS-251-18	18 cm (7")

Reynolds Dissecting Scissor
Curved



Catalog #	Size
SS-251-15 SC	15.5 cm (6")
SS-251-18 SC	18 cm (7")

Jameson Dissecting Scissor
Curved



Catalog #	Size
SS-253-15	15.5 cm (6")


Jameson Dissecting Scissor
Curved



SUPER Edge™

Catalog #	Size
SS-253-15 SC	15.5 cm (6")

Jameson-Werber Dissecting Scissor
Curved




Catalog #	Size
SS-257-13	13 cm (5")



Dissecting Scissors, Vascular Scissors
Präparierscheren, Gefäßscheren
Tijeras de disección, Tijeras vasculares
Ciseaux à dissection, Ciseaux vasculaires
Forbici per dissezione, Forbici vascolari


Potts Dissecting Scissor
Curved

Catalog #	Size
SS-259-18	18 cm (7")




DeBakey Arteriotomy Scissor
Curved

Catalog #	Size
SS-261-17	17.5 cm (7")




Vascular Scissor
Straight

Catalog #	Size
SS-264-19	19 cm (7 1/2")




Vascular Scissor
Curved

Catalog #	Size
SS-265-19	19 cm (7 1/2")



Satinsky Thoracic Scissor
Curved

Catalog #	Size
SS-267-24	25.5 cm (10")






Vascular Scissors
Gefäßscheren
Tijeras vasculares
Ciseaux vasculaires
Forbici vascolari

SS

SCISSORS

Diethrich-Salyer Cleft Palate Dissecting Scissor
Angled

Catalog #	Size
SS-269-12	12 cm (4 ¾")



Diethrich Vascular Scissor
Angled 125°, With Guide

Catalog #	Size
SS-275-00	17.5 cm (7")



DeBakey Vascular Scissor
Angled 45°

Catalog #	Size
SS-277-45	16.5 cm (6 ½")



DeBakey Vascular Scissor
Angled 60°

Catalog #	Size
SS-277-60	16cm (6 ¼")



Potts-Smith Vascular Scissor
Angled 60°, With Probe Tip

Catalog #	Size
SS-279-19	18 cm (7")






Vascular Scissors
Gefäßscheren
Tijeras vasculares
Ciseaux vasculaires
Forbici vascolari


Potts-Smith Vascular Scissor
Angled 25°

Catalog #	Size
SS-281-25	18 cm (7")




Potts-Smith Vascular Scissor
Angled 25°

Catalog #	Size
SS-281-25 SC	18 cm (7")





Potts-Smith Vascular Scissor
Angled 45°

Catalog #	Size
SS-281-45	18.5 cm (7 1/4")





Potts-Smith Vascular Scissor
Angled 45°

Catalog #	Size
SS-281-45 SC	18.5 cm (7 1/4")



Potts-Smith Vascular Scissor
Angled 60°

Catalog #	Size
SS-281-60	18 cm (7")






Vascular Scissors
Gefäßscheren
Tijeras vasculares
Ciseaux vasculaires
Forbici vascolari

SS

SCISSORS

Potts-Smith Vascular Scissor
Angled 60°

Catalog #	Size
SS-281-60 SC	18 cm (7")



Potts-Smith Vascular Scissor
Angled on Flat

Catalog #	Size
SS-283-19	18.5 cm (7 1/4")





Vascular Scissors
Gefäßscheren
Tijeras vasculares
Ciseaux vasculaires
Forbici vascolari

Diethrich-Potts Vascular Scissor
Standard Blade



Catalog #	Size, Profile
SS-285-25	18 cm (7"), Angled 25°
SS-285-45	18 cm (7"), Angled 45°
SS-285-60	18 cm (7"), Angled 60°
SS-285-90	18 cm (7"), Angled 90°
SS-287-00	18 cm (7"), Angled 125°




SS-285-25, Angled 25°




SS-285-45, Angled 45°



SS-285-60, Angled 60°



SS-285-90, Angled 90°



SS-287-00, Angled 125°



Vascular Scissors
Gefäßscheren
Tijeras vasculares
Ciseaux vasculaires
Forbici vascolari

SS

SCISSORS

Diethrich-Potts Vascular Scissor
Delicate Blade




Catalog #	Size, Profile
SS-289-25	18 cm (7"), Angled 25°
SS-289-45	18 cm (7"), Angled 45°
SS-289-60	18 cm (7"), Angled 60°
SS-289-90	18 cm (7"), Angled 90°
SS-291-00	18 cm (7"), Angled 125°




SS-289-25, Angled 25°




SS-289-45, Angled 45°



SS-289-60, Angled 60°



SS-289-90, Angled 90°



SS-291-00, Angled 125°



Vascular Scissors
Gefäßscheren
Tijeras vasculares
Ciseaux vasculaires
Forbici vascolari

Diethrich-Potts Vascular Scissor
Ultra Delicate Blade



Catalog #	Size, Profile
SS-293-25	18 cm (7"), Angled 25°
SS-293-45	18 cm (7"), Angled 45°
SS-293-60	18 cm (7"), Angled 60°
SS-293-90	18 cm (7"), Angled 90°
SS-295-00	18 cm (7"), Angled 125°




SS-293-25, Angled 25°




SS-293-45, Angled 45°



SS-293-60, Angled 60°



SS-293-90, Angled 90°



SS-295-00, Angled 125°



Vascular Scissors
Gefäßscheren
Tijeras vasculares
Ciseaux vasculaires
Forbici vascolari

SS

SCISSORS

Potts-De Martel Cardiovascular Scissor
Angled 60°

Catalog #	Size
SS-299-22	22 cm (8 ¾")




DeBakey Offset Scissor
S Shaped

Catalog #	Size
SS-301-16	15.5 cm (6")




Favaloro Offset Scissor
Angled

Catalog #	Size
SS-303-15	15 cm (6")




Vascular Scissor
One Blade Probe

Catalog #	Size
SS-306-12	12 cm (4 ¾")



Vascular Scissor
One Blade Probe



Catalog #	Size
SS-308-14	14.5 cm (5 ¾")






Vascular Scissors
Gefäßscheren
Tijeras vasculares
Ciseaux vasculaires
Forbici vascolari

Strully Vascular Scissor
Curved




Catalog #	Size
SS-311-22	22 cm (8 3/4")

Vascular Scissor
Straight




Catalog #	Size
SS-312-22	22.5 cm (8 3/4")

Barre Vascular Scissor
Straight





Catalog #	Size
SS-314-28	28 cm (11")

Potts-Smith Vascular Scissor
Curved



Catalog #	Size
SS-315-19	19 cm (7 1/2")

Cooley Vascular Scissor
Curved



Catalog #	Size
SS-317-19	19 cm (7 1/2")




Vascular Scissors, Dissecting Scissors
Gefäßscheren
Tijeras vasculares, Tijeras de disección
Ciseaux vasculaires, Ciseaux à dissection
Forbici vascolari, Forbici per dissezione

SS

SCISSORS

Klinkenberg-Loth Vascular Scissor
S Shaped

Catalog #	Size
SS-319-23	23 cm (9")



Olivecrona Trigeminal Scissor
Toothed

Catalog #	Size, Working Le
SS-323-23	23 cm (9"), 120




Strully Neurological Scissor
Curved

Catalog #	Size
SS-325-23	22.5 cm (8 3/4")



Schmieden-Taylor Dissecting Scissor
One Blade Probed

Catalog #	Size
SS-327-17	16.5 cm (6 1/2")





Dissecting Scissors
Präparierscheren
Tijeras de disección
Ciseaux à dissection
Forbici per dissezione


Toennis-Adson Dissecting Scissor
Curved

Catalog #	Size
SS-329-17	17.5 cm (7")




Toennis Dissecting Scissor
Straight, Very Delicate

Catalog #	Size
SS-332-18	18 cm (7")



Toennis Dissecting Scissor
Straight, Delicate

Catalog #	Size
SS-334-18	18 cm (7")




Toennis Dissecting Scissor
Curved, Delicate

Catalog #	Size
SS-335-18	18 cm (7")



Killian Bone Shear

Catalog #	Size
SS-341-21	21 cm (8 1/4")






Tonsil and Nasal Scissors
Tonsillen und Nasenscheren
Tijeras para amígdalas, Tijeras nasales
Ciseaux pour les amygdales, Ciseaux nasaux
Forbici per tonsille, Forbici nasali

SS

SCISSORS

Muehling Bone Shear

Catalog #	Size
SS-343-21	21 cm (8 1/4")



Caplan Septum Scissor
One Toothed Cutting Edge

Catalog #	Size
SS-345-20	20 cm (8")




Prince Tonsil Scissor

Catalog #	Size
SS-347-18	18 cm (7")




Dean Tonsil Scissor
Angled Upwards

Catalog #	Size
SS-351-17	17.5 cm (7")



Dean Tonsil Scissor
Toothed, Angled Upwards


Catalog #	Size
SS-353-17	17.5 cm (7")





Tonsil and Nasal Scissors
Tonsillen und Nasenscheren
Tijeras para amígdalas, Tijeras nasales
Ciseaux pour les amygdales, Ciseaux nasaux
Forbici per tonsille, Forbici nasali

Good Tonsil Scissor



Catalog #	Size
SS-355-19	19.5 cm (7 3/4")

Fomon Nasal Scissor



Catalog #	Size
SS-357-15	15 cm (6")



Fomon Nasal Scissor



SUPER Edge™

Catalog #	Size
SS-357-15 SC	15 cm (6")

Cottle Nasal Scissor




Catalog #	Size
SS-359-16	16 cm (6 1/4")

Cottle Nasal Scissor

SUPER Edge™

Catalog #	Size
SS-359-16 SC	16 cm (6 1/4")



Nasal Scissors, Episiotomy Scissors
Nasenscheren, Episiotomiescheren
Tijeras nasals, Tijeras para episiotomía
Ciseaux nasaux, Ciseaux à épisiotomie
Forbici nasali, Forbici per episiotomia

SS
SCISSORS
Heymann Nasal Scissor

Catalog #	Size
SS-361-18	18.5 cm (7 ¼")


Heymann Nasal Scissor

Catalog #	Size
SS-361-18 SC	18.5 cm (7 ¼")


Heymann Nasal Scissor
One Toothed Cutting Edge

Catalog #	Size
SS-363-18	18.5 cm (7 ¼")


Converse (Graefe) Nasal Scissor
Angled, Blunt/Blunt

Catalog #	Size
SS-367-11	10.5 cm (4 ¼")


Converse Nasal Scissor
Angled, Sharp/Sharp

Catalog #	Size
SS-369-11	10.5 cm (4 ¼")






Nasal Scissors, Episiotomy Scissors, Umbilical Scissors
Nasenscheren, Episiotomiescheren, Nabelschnurscheren
Tijeras nasals, Tijeras para episiotomía, Tijeras umbilicales
Ciseaux nasaux, Ciseaux à épisiotomie, Ciseaux pour le cordon ombilical
Forbici nasali, Forbici per episiotomia, Forbici ombelicali

Waldmann Nasal Scissor
One Toothed Cutting Edge

Catalog #	Size
SS-373-18	18 cm (7")



Braun-Stadler Episiotomy Scissor

Catalog #	Size
SS-381-14	14.5 cm (5 3/4")



Braun-Stadler Episiotomy Scissor

Catalog #	Size
SS-383-20	20 cm (8")
SS-383-22	22 cm (8 3/4")




Modell USA Umbilical Cord Scissor

Catalog #	Size
SS-385-10	10.5 cm (4 1/4")



Schumacher Umbilical Cord Scissor

Catalog #	Size
SS-387-16	16 cm (6 1/4")






Umbilical Scissors, Gum Scissors
Nabelschnurscheren, Zahnfleischscheren
Tijeras umbilicales, Tijeras para encías
Ciseaux pour le cordon ombilical, Ciseaux à gencives
Forbici ombelicali, Forbici per gengiva

SS

SCISSORS


Busch Umbilical Cord Scissor

Catalog #	Size
SS-389-16	16 cm (6 ¼")




Busch Umbilical Cord Scissor
One Toothed Cutting Edge

Catalog #	Size
SS-391-16	16 cm (6 ¼")



Goldman-Fox Gum Scissor
Straight


Catalog #	Size
SS-400-13	13 cm (5")



Goldman-Fox Gum Scissor
Straight


SUPER *Edge*™

Catalog #	Size
SS-400-13 SC	13 cm (5")



Goldman-Fox Gum Scissor
Curved

Catalog #	Size
SS-401-13	13 cm (5")






Gum Scissors, Operating Scissors
Zahnfleischscheren, Chirurgische Scheren
Tijeras para encías, Tijeras quirúrgicas
Ciseaux à gencives, Ciseaux chirurgicaux
Forbici per gengiva, Forbici per chirurgia

Goldman-Fox Gum Scissor
Curved

SUPER Edge™

Catalog #	Size
SS-401-13 SC	13 cm (5")



Wagner Operating Scissor
Straight

Catalog #	Size
SS-402-12	12 cm (4 3/4")



Wagner Operating Scissor
Straight

SUPER Edge™

Catalog #	Size
SS-402-12 SC	12 cm (4 3/4")



Wagner Operating Scissor
Curved

Catalog #	Size
SS-403-12	12 cm (4 3/4")



Wagner Operating Scissor
Curved

SUPER Edge™

Catalog #	Size
SS-403-12 SC	12 cm (4 3/4")




Operating Scissors, Gum Scissors
Chirurgische Scheren, Zahnfleischscheren
Tijeras quirúrgicas, Tijeras para encías
Ciseaux chirurgicaux, Ciseaux à gencives
Forbici per chirurgia, Forbici per gengiva

SS
SCISSORS


Fine Pattern Operating Scissor
Straight, Sharp/Blunt

Catalog #	Size
SS-406-12	12 cm (4 ³ / ₄ ")




Fine Pattern Operating Scissor
Curved, Sharp/Blunt

Catalog #	Size
SS-407-12	12 cm (4 ³ / ₄ ")




Fine Pattern Operating Scissor
Straight, Sharp/Sharp

Catalog #	Size
SS-408-12	12 cm (4 ³ / ₄ ")




Fine Pattern Operating Scissor
Curved, Sharp/Sharp

Catalog #	Size
SS-409-12	12 cm (4 ³ / ₄ ")



La Grange Gum Scissor
S Shaped

Catalog #	Size
SS-411-11	11.5 cm (4 ¹ / ₂ ")





Gum Scissors, Operating Scissors
Zahnfleischscheren, Chirurgische Scheren
Tijeras para encías, Tijeras quirúrgicas
Ciseaux à gencives, Ciseaux chirurgicaux
Forbici per gengiva, Forbici per chirurgia


Delicate Gum Scissor
Straight

Catalog #	Size
SS-414-11	11.5 cm (4 ½")




Delicate Gum Scissor
Curved

Catalog #	Size
SS-415-11	11.5 cm (4 ½")




Locklin Gum Scissor
Angled, One Toothed Cutting Edge

Catalog #	Size
SS-417-16	16 cm (6 ¼")




Locklin Gum Scissor
Angled, S Shaped, One Toothed Cutting Edge

Catalog #	Size
SS-419-16	16 cm (6 ¼")



Very Delicate Operating Scissor
Straight

Catalog #	Size
SS-422-09	9 cm (3 ½")






Operating Scissors, Iris Scissors
Chirurgische Scheren, Irisscheren
Tijeras quirúrgicas, Tijeras para iridectomía
Ciseaux chirurgicaux, Ciseaux à iridectomie
Forbici per chirurgia, Forbici per iridectomia


Very Delicate Operating Scissor
Curved

Catalog #	Size
SS-423-09	9 cm (3 ½")




Iris Scissor
Straight

Catalog #	Size
SS-430-09	9 cm (3 ½")
SS-430-10	10 cm (4")




Iris Scissor
Curved

Catalog #	Size
SS-431-09	9 cm (3 ½")
SS-431-10	10 cm (4")




Iris Scissor
With Extra Flat Blades

Catalog #	Size
SS-432-10	10 cm (4")



Bonn Iris Scissor
Straight

Catalog #	Size
SS-434-09	9 cm (3 ½")






Iris Scissors
Irisscheren
Tijeras para iridectomía
Ciseaux à iridectomie
Forbici per iridectomia


Bonn Iris Scissor
Curved

Catalog #	Size
SS-435-09	9 cm (3 ½")



Iris Scissor
Straight

Catalog #	Size
SS-440-11	11.5 cm (4 ½")



Iris Scissor
Straight


SUPER Edge™

Catalog #	Size
SS-440-11 SC	11.5 cm (4 ½")



Iris Scissor
Curved

Catalog #	Size
SS-441-11	11.5 cm (4 ½")



Iris Scissor
Curved

Catalog #	Size
SS-441-11 SC	11.5 cm (4 ½")






Iris Scissors
Irisscheren
Tijeras para iridectomía
Ciseaux à iridectomie
Forbici per iridectomia

SS

SCISSORS


Iris Scissor
Curved to Side

Catalog #	Size
SS-443-11	11.5 cm (4 ½")



Iris Scissor
Angled

Catalog #	Size
SS-445-11	11.5 cm (4 ½")




Iris Scissor
Straight

Catalog #	Size
SS-446-12	12.5 cm (5")




Iris Scissor
Curved

Catalog #	Size
SS-447-12	12.5 cm (5")



Stevens Tenotomy Scissor
Straight, Sharp/Sharp

Catalog #	Size
SS-450-10	10 cm (4")
SS-450-11	11 cm (4 ½")






Iris Scissors, Tendon Scissors
Irisscheren, Sehnenscheren
Tijeras para iridectomía, Tijeras para tendones
Ciseaux à iridectomie, Ciseaux à ténotomie
Forbici per iridectomia, Forbici per tenotomia

Stevens Tenotomy Scissor
Straight, Sharp/Sharp


SUPER Edge™

Catalog #	Size
SS-450-11 SC	11 cm (4 ½")



Stevens Tenotomy Scissor
Curved, Sharp/Sharp


Catalog #	Size
SS-451-10	10 cm (4")
SS-451-11	11 cm (4 ½")



Stevens Tenotomy Scissor
Curved, Sharp/Sharp


SUPER Edge™

Catalog #	Size
SS-451-11 SC	11 cm (4 ½")



Stevens Tenotomy Scissor
Straight, Blunt/Blunt


Catalog #	Size
SS-452-10	10 cm (4")
SS-452-11	11 cm (4 ½")
SS-452-12	12.5 cm (5")



Stevens Tenotomy Scissor
Straight, Blunt/Blunt

SUPER Edge™

Catalog #	Size
SS-452-11 SC	11 cm (4 ½")






Iris Scissors, Strabismus Scissors
Irisscheren, Strabismusscheren
Tijeras para iridectomia, Tijeras para estrabismo
Ciseaux à iridectomie, Ciseaux à strabisme
Forbici per iridectomia, Forbici per strabismo

SS

SCISSORS

Stevens Tenotomy Scissor
Curved, Blunt/Blunt


Catalog #	Size
SS-453-10	10 cm (4")
SS-453-11	11 cm (4 ½")
SS-453-12	12.5 cm (5")



Stevens Tenotomy Scissor
Curved, Blunt/Blunt


SUPER Edge™

Catalog #	Size
SS-453-11 SC	11 cm (4 ½")



Strabismus Scissor
Straight


Catalog #	Size
SS-456-11	11.5 cm (4 ½")



Strabismus Scissor
Straight


SUPER Edge™

Catalog #	Size
SS-456-11 SC	11.5 cm (4 ½")



Strabismus Scissor
Curved

Catalog #	Size
SS-457-11	11.5 cm (4 ½")






Strabismus Scissors, Dissecting Scissors
Strabismusscheren, Präparierscheren
Tijeras para estrabismo, Tijeras de disección
Ciseaux à strabisme, Ciseaux à dissection
Forbici per strabismo, Forbici per dissezione

Strabismus Scissor
Curved


SUPER Edge™

Catalog #	Size
SS-457-11 SC	11.5 cm (4 1/2")



Littler Dissecting Scissor
Curved, Tip With Eye For Suture

Catalog #	Size
SS-461-11	11.5 cm (4 1/2")



Cottle Bulldog Scissor
Curved

Catalog #	Size
SS-463-11	11.5 cm (4 1/2")




Kilner Delicate Dissecting Scissor
Straight

Catalog #	Size
SS-466-11	11.5 cm (4 1/2")



Enucleation Scissor
Curved


Catalog #	Size
SS-469-11	11.5 cm (4 1/2")




Dissecting Scissors, Face-lift Scissors
 Präparierscheren, Face-lift-Scheren
 Tijeras de disección, Tijeras para Face-Lift
 Ciseaux à dissection, Ciseaux pour Face-lift
 Forbici per dissezione, Forbici per Face-lift

SS
SCISSORS
Chadwick Delicate Scissor
Curved


Catalog #	Size
SS-471-11	11.5 cm (4 ½")


Quimby Delicate Scissor
Curved


Catalog #	Size
SS-473-13	13 cm (5")


Nerve Dissecting Scissor
Straight


Catalog #	Size
SS-476-11	11 cm (4 ½")
SS-476-13	13 cm (5")


Nerve Dissecting Scissor
Curved

Catalog #	Size
SS-477-11	11 cm (4 ½")
SS-477-13	13 cm (5")


Delicate Scissor
Angled

Catalog #	Size
SS-481-11	11.5 cm (4 ½")






Face-lift Scissors
Face-lift-Scheren
Tijeras para Face-Lift
Ciseaux pour Face-lift
Forbici per Face-lift


Ribbon Salyer Scissor
Straight, Flat Shanks

Catalog #	Size
SS-484-09	9.5 cm (3 3/4")




Ribbon Iris Scissor
Straight, Sharp, Flat Shanks

Catalog #	Size
SS-486-09	9.5 cm (3 3/4")




Ribbon Iris Scissor
Curved, Sharp, Flat Shanks

Catalog #	Size
SS-487-09	9.5 cm (3 3/4")




Ribbon Iris Scissor
Straight, Blunt, Flat Shanks

Catalog #	Size
SS-488-10	10 cm (4")



Ribbon Iris Scissor
Curved, Blunt, Flat Shanks


Catalog #	Size
SS-489-10	10 cm (4")




Face-lift Scissors
Face-lift-Scheren
Tijeras para Face-Lift
Ciseaux pour Face-lift
Forbici per Face-lift

SS
SCISSORS
Ribbon Strabismus Scissor
Straight, Flat Shanks


Catalog #	Size
SS-492-10	10 cm (4")


Ribbon Strabismus Scissor
Curved, Flat Shanks


Catalog #	Size
SS-493-10	10 cm (4")


Ribbon Stevens Tenotomy Scissor
Straight, Flat Shanks, Sharp Tips


Catalog #	Size
SS-496-10	10 cm (4")


Ribbon Stevens Tenotomy Scissor
Curved, Flat Shanks, Sharp Tips

Catalog #	Size
SS-497-10	10 cm (4")


Ribbon Stevens Tenotomy Scissor
Straight, Flat Shanks, Blunt Tips

Catalog #	Size
SS-498-10	10 cm (4")






Face-lift Scissors, Enterotomy Scissors
Face-lift-Scheren, Darmscheren
Tijeras para Face-Lift, Tijeras para intestino
Ciseaux pour Face-lift, Ciseaux-entérotomes
Forbici per Face-lift, Forbici per enterotomia


Ribbon Stevens Tenotomy Scissor
Curved, Flat Shanks, Blunt Tips

Catalog #	Size
SS-499-10	10 cm (4")




Ribbon Gradle Scissor
Straight, Flat Shanks

Catalog #	Size
SS-500-10	10 cm (4")




Ribbon Gradle Scissor
Curved, Flat Shanks

Catalog #	Size
SS-501-10	10 cm (4")




Ribbon Kaye Face-Lift Scissor
Curved, Flat Shanks, Toothed

Catalog #	Size
SS-507-11	11 cm (4 1/2")



Enterotomy Scissor
One Probe Tip

Catalog #	Size
SS-520-21	21.5 cm (8 1/2")






Rectal Scissors
Rektumschere
Tijera para recto
Ciseaux pour rectum
Forbici per il retto

SS


SCISSORS

Enterotomy Scissor
One Probe Tip



Catalog #	Size
SS-522-21	21.5 cm (8 ½")

Rectal Scissor
Curved, Heavy Pattern




Catalog #	Size
SS-525-32	32.5 cm (12 ¾")

Bowel Dissecting Scissor
Angled 60°




Catalog #	Size
SS-527-28	28 cm (11")

Lloyd-Davies (Goligher) Rectal Scissor
Curved



Catalog #	Size
SS-529-27	27 cm (10 ¾")

Stelzner Rectal Scissor
Straight




Catalog #	Size
SS-530-31	31 cm (12 ¼")



Rectal Scissors, Stitch Scissors
Rektumschere, Fadenscheren
Tijera para recto, Tijeras para hilos
Ciseaux pour rectum, Ciseaux à ligatures
Forbici per il retto, Forbici per legature


Stelzner Rectal Scissor
Curved

Catalog #	Size
SS-531-31	31 cm (12 1/4")




Mueller Rectal Scissor
S Shaped

Catalog #	Size
SS-533-32	32.5 cm (12 3/4")




Northbent Stitch Scissor

Catalog #	Size
SS-540-13	13 cm (5")




Spencer Stitch Scissor

Catalog #	Size
SS-542-09	9 cm (3 1/2")
SS-542-13	13 cm (5")



Buck Stitch Scissor


Catalog #	Size
SS-544-14	14.5 cm (5 3/4")
SS-544-18	18 cm (7")



Stitch Scissors, Wire Cutting Scissors, Crown Scissors
Fadenscheren, Drahtscheren, Kronenscheren
Tijeras para hilos, Tijeras para cortar hilos de acero, Tijeras para coronas
Ciseaux à ligatures, Ciseaux coupe-fil, Ciseaux à couronnes
Forbici per legature, Forbici per tagliare fili metallici, Forbici per corone

Littauer Stitch Scissor



Catalog #	Size
SS-546-14	14 cm (5 1/2")

Eiselberg Stitch Scissor



Catalog #	Size
SS-548-12	12 cm (4 3/4")

Universal Wire Cutting Scissor
Angled, One Toothed Cutting Edge


Catalog #	Size
SS-551-12	12 cm (4 3/4")

Harvey Wire Cutting Scissor
Straight, One Toothed Cutting Edge


Catalog #	Size
SS-554-12	12.5 cm (5")

Beebee Crown Scissor
Straight



Catalog #	Size
SS-556-10	10.5 cm (4 1/4")



Crown Scissors, Nail Scissors
Kronenscheren, Nagelscheren
Tijeras para coronas, Tijeras para uñas
Ciseaux à couronnes, Ciseaux à ongles
Forbici per corone, Forbici per unghie


Beebee Crown Scissor
Curved

Catalog #	Size
SS-557-10	10.5 cm (4 1/4")




Beebee Crown Scissor
Straight, One Toothed Cutting Edge

Catalog #	Size
SS-558-12	12 cm (4 3/4")



Beebee Crown Scissor
Curved, One Toothed Cutting Edge

Catalog #	Size
SS-559-12	12 cm (4 3/4")




Nail Scissor

Catalog #	Size
SS-571-09	9 cm (3 1/2")



Nail Scissor
Delicate Pattern

Catalog #	Size
SS-573-09	9 cm (3 1/2")






Nail Scissors
Nagelscheren
Tijeras para uñas
Ciseaux à ongles
Forbici per unghie

SS

SCISSORS


Nail Scissor

Catalog #	Size
SS-575-09	8.5 cm (3 1/4")



**Nail Scissor
Straight**

Catalog #	Size
SS-576-09	9 cm (3 1/2")




**Nail Scissor
Curved**

Catalog #	Size
SS-577-09	9 cm (3 1/2")
SS-577-10	10 cm (4")
SS-577-11	11 cm (4 1/2")




**Nail Scissor
Straight**

Catalog #	Size
SS-578-09	9.5 cm (3 3/4")



**Nail Scissor
Curved**

Catalog #	Size
SS-579-09	9.5 cm (3 3/4")
SS-579-10	10.5 cm (4 1/4")






Bandage Scissors
Verbandscheren
Tijeras para vendajes
Ciseaux pour pansements
Forbici per bendaggi


Smith (Mod. USA) Bandage Scissor

Catalog #	Size
SS-582-18	18 cm (7")
SS-582-20	20 cm (8")
SS-583-23	23 cm (9")



**Bandage Scissor
Straight**

Catalog #	Size
SS-584-15	15 cm (6")



**Lorenz Bandage Scissor
Heavy Pattern**

Catalog #	Size
SS-585-23	23 cm (9")



**Knowles Bandage Scissor
Angled**

Catalog #	Size
SS-587-14	14 cm (5 1/2")



**Knowles Bandage Scissor (For Fingers)
Straight**

Catalog #	Size
SS-588-14	14cm (5 1/2")





Bandage Scissors
Verbandscheren
Tijeras para vendajes
Ciseaux pour pansements
Forbici per bendaggi

SS

SCISSORS


Lister Bandage Scissor

Catalog #	Size
SS-591-14	14.5 cm (5 ¾")
SS-591-18	18.5 cm (7 ¼")
SS-591-20	20 cm (8")



Lister Bandage Scissor
One Large Ring

Catalog #	Size
SS-593-18	18.5 cm (7 ¼")
SS-593-21	21 cm (8 ¼")



Lister (Eccentric) Bandage Scissor
One Toothed Cutting Edge

Catalog #	Size
SS-595-16	16 cm (6 ¼")
SS-595-19	19 cm (7 ½")



Bergmann Bandage Scissor
Heavy Pattern


Catalog #	Size
SS-597-23	23.5 cm (9 ¼")





Bandage Scissors, Plaster Shears
Verbandscheren, Gipsscheren
Tijeras para vendajes, Tijeras para yeso
Ciseaux pour pansements, Ciseaux pour plâtre
Forbici per bendaggi, Forbici per gesso

Esmarch Bandage Scissor
Heavy Pattern, Short Blades



Catalog #	Size
SS-599-20	20 cm (8")
SS-599-23	23 cm (9")

Universal Bandage Scissor
Plastic Handle



Catalog #	Size	Handle Colour
SS-601-14	14.5 cm (5 ¾")	Black
SS-601-16	16 cm (6 ¼")	Black
SS-601-18	18 cm (7 ½")	Black
SS-603-14	14.5 cm (5 ¾")	Blue
SS-603-16	16 cm (6 ¼")	Blue
SS-603-18	18 cm (7 ½")	Blue
SS-605-14	14.5 cm (5 ¾")	Red
SS-605-16	16 cm (6 ¼")	Red
SS-605-18	18 cm (7 ½")	Red
SS-607-14	14.5 cm (5 ¾")	Green
SS-607-16	16 cm (6 ¼")	Green
SS-607-18	18 cm (7 ½")	Green

Seutin Plaster of Paris Shear



Catalog #	Size
SS-609-23	23 cm (9")



Plaster Shears
Gipsscheren
Tijeras para yeso
Ciseaux pour plâtre
Forbici per gesso

SS

SCISSORS


Plaster of Paris Shear
One Toothed Cutting Edge

Catalog #	Size
SS-611-20	20 cm (8")



Bruns Plaster of Paris Shear
One Toothed Cutting Edge

Catalog #	Size
SS-613-24	24 cm (9 1/2")





Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

Noyes Iris Scissor
Straight, Sharp/Sharp

Catalog #	Size
SS-630-12	12.5 cm (5")



Noyes Iris Scissor
Curved, Sharp/Sharp

Catalog #	Size
SS-631-12	12.5 cm (5")



Noyes Iris Scissor
Straight, Sharp/Blunt

Catalog #	Size
SS-632-12	12.5 cm (5")



Noyes Iris Scissor
Curved, Sharp/Blunt

Catalog #	Size
SS-633-12	12.5 cm (5")



Noyes Iris Scissor
Straight, Blunt/Blunt

Catalog #	Size
SS-634-12	12.5 cm (5")





Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

SS

SCISSORS

Noyes Iris Scissor
Curved, Blunt/Blunt



Catalog #	Size
SS-635-12	12.5 cm (5")

Micro Scissor
Straight



Catalog #	Size
SS-636-12	12 cm (4 3/4")

Micro Scissor
Curved



Catalog #	Size
SS-637-12	12 cm (4 3/4")

Wecker Iris Scissor
Angled, Sharp/Sharp



Catalog #	Size
SS-641-11	11 cm (4 1/2")

Wecker Iris Scissor
Angled, Sharp/Blunt



Catalog #	Size
SS-643-11	11 cm (4 1/2")



Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

Wecker Iris Scissor
Angled, Blunt/Blunt



Catalog #	Size
SS-645-11	11 cm (4 1/2")

Barraquer Iris Scissor
Angled



Catalog #	Size, Blade Size
SS-647-05	5.5 cm (2 1/4"), 5 mm
SS-647-07	5.5 cm (2 1/4"), 7 mm
SS-647-10	5.5 cm (2 1/4"), 10 mm

Castroviejo Corneoscleral Scissor
Small Blades



Catalog #	Size, Profile
SS-651-01	10.5 cm (4"), Right
SS-651-02	10.5 cm (4"), Left

Castroviejo Corneoscleral Scissor
Medium Blades



Catalog #	Size
SS-653-01	11 cm (4 1/2"), Right
SS-653-02	11 cm (4 1/2"), Left



Troutman-Castroviejo Corneoscleral Scissor

Small Blades, With Lock

Catalog #	Size, Profile
SS-655-01	10.5 cm (4"), Right
SS-655-02	10.5 cm (4"), Left



Troutman-Castroviejo Corneoscleral Scissor

Medium Blades, With Lock

Catalog #	Size, Profile
SS-657-01	11 cm (4 1/2"), Right
SS-657-02	11 cm (4 1/2"), Left



Troutman-Katzin Corneal Transplant Scissor

Strongly Curved, Small Blades With Lock

Catalog #	Size, Profile
SS-659-01	10.5 cm (4"), Right
SS-659-02	10.5 cm (4"), Left



Katzin Corneal Transplant Scissor

Strongly Curved, Medium Blades With Lock

Catalog #	Size, Profile
SS-661-01	11 cm (4 1/2"), Right
SS-661-02	11 cm (4 1/2"), Left





Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

Castroviejo Universal Cornel Scissor
Curved, Blunt Tips, Small Blades

Catalog #	Size, Profile
SS-663-11	11 cm (4 1/2"), Small Blades
SS-665-11	11 cm (4 1/2"), Medium Blades
SS-667-11	11 cm (4 1/2"), Large Blades



Castroviejo Universal Cornel Scissor
Strongly Curved, Blunt Tips

Catalog #	Size
SS-669-11	11.5 cm (4 1/2")



Castroviejo Keratoplasty Scissor
Angled, Blunt Tips, Medium Blades

Catalog #	Size
SS-671-11	11.5 cm (4 1/2")



McPherson-Westcott Conjunctival Scissor
Curved, Blunt Tips, Small Blades

Catalog #	Size
SS-673-11	11 cm (4 1/2")



Westcott Tenotomy Scissor
Curved, Blunt Tips, Medium Blades

Catalog #	Size
SS-675-11	11 cm (4 1/2")





Westcott Tenotomy Scissor

Curved, Blunt Tips, Standard Blades

Catalog #	Size, Profile
SS-677-01	11.5 cm (4 ½"), Right
SS-677-02	11.5 cm (4 ½"), Left



Tibolt Punctal Scissor

One Conically Shaped Blade with Micro Grooves & Textured Surface

Catalog #	Size
SS-679-11	11.5 cm (4 ½")



McPherson-Westcott Stitch Scissor

Curved, Very Sharp Pointed Tips, Small Blades

Catalog #	Size
SS-681-11	11 cm (4 ½")



Westcott Stitch Scissor

Curved, Very Sharp Pointed Tips, Standard Blades

Catalog #	Size
SS-683-11	11.5 cm (4 ½")



Jaffe Stitch Scissor

Curved, Very Sharp Pointed Tips, Medium Blades

Catalog #	Size
SS-685-11	11 cm (4 ½")





Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

Iris Micro Scissor

Straight, Sharp Tips, Small Blades

Catalog #	Size
SS-688-11	11 cm (4 1/2")



Iris Micro Scissor

Curved, Sharp Tips, Small Blades

Catalog #	Size
SS-689-11	11 cm (4 1/2")



Fine Micro Scissor

Straight, Sharp Tips, Medium Blades

Catalog #	Size
SS-690-11	11.5 cm (4 1/2")



Fine Stitch Scissor

Curved, Sharp Tips, Medium Blades

Catalog #	Size
SS-691-11	11.5 cm (4 1/2")



Vannas Capsulotomy Scissor

Straight, Sharp Tips

Catalog #	Size, Blade Size
SS-694-08	8 cm (3 1/4"), 5 mm



Vannas Capsulotomy Scissor

Curved, Sharp Tips

Catalog #	Size, Blade Size
SS-695-08	8 cm (3 1/4"), 5 mm


Vannas Capsulotomy Scissor

Angled Forward, Sharp Tips

Catalog #	Size, Blade Size
SS-696-08	8 cm (3 1/4"), 5 mm


Vannas Capsulotomy Scissor

Angled to Side, Sharp Tips

Catalog #	Size, Blade Size
SS-697-08	8 cm (3 1/4"), 5 mm


Gills-Vannas Capsulotomy Scissor

Straight, Sharp Tips

Catalog #	Size, Blade Size
SS-698-08	8 cm (3 1/4"), 7 mm


Gills-Vannas Capsulotomy Scissor

Curved, Sharp Tips

Catalog #	Size, Blade Size
SS-699-08	8 cm (3 1/4"), 7 mm





Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

Gills-Vannas Capsulotomy Scissor
Angled Forward, Sharp Tips



Catalog #	Size, Blade Size
SS-701-08	8 cm (3 1/4"), 7 mm

McPherson-Vannas Iris Scissor
Curved, Sharp Tips



Catalog #	Size
SS-703-08	8 cm (3 1/4")

Gills-Vannas Micro Scissor
Straight, Sharp Tips, Extra Thin



Catalog #	Size, Blade Size
SS-706-10	10.5 cm (4"), 7 mm

Gills-Vannas Micro Scissor
Curved, Sharp Tips, Extra Thin



Catalog #	Size, Blade Size
SS-707-10	10.5 cm (4"), 7 mm

Gills-Vannas Micro Scissor
Angled Forward, Sharp Tips, Extra Thin



Catalog #	Size, Blade Size
SS-709-10	10.5 cm (4"), 7 mm



Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

SS

SCISSORS

Stern-Gills Micro Scissor

Straight, Sharp Tips, Extra Thin



Catalog #	Size, Blade Size
SS-710-11	11 cm (4 1/2"), 10 mm

Stern-Gills Micro Scissor

Curved, Sharp Tips, Extra Thin



Catalog #	Size, Blade Size
SS-711-11	11 cm (4 1/2"), 10 mm

Stern-Gills Micro Scissor

Angled Forward, Sharp Tips, Extra Thin



Catalog #	Size, Blade Size
SS-713-11	11 cm (4 1/2"), 10 mm

Gills Capsulotomy Scissor


Angled Blades With Curved Tips, Sharp Points



Catalog #	Size, Blade Size
SS-715-11	10.5 cm (4"), 10 mm

Vannas Micro Scissor

Straight



Catalog #	Size
SS-726-12	12 cm (4 3/4")



Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

Vannas Micro Scissor
Curved

Catalog #	Size
SS-727-12	12 cm (4 3/4")



Vannas Micro Scissor
Angled on Flat

Catalog #	Size
SS-729-12	12 cm (4 3/4")




Biemer Vein Scissor

Catalog #	Size
SS-730-13	13 cm (5")




Micro Scissor
Straight

Catalog #	Size
SS-732-16	16 cm (6 1/4")



Micro Scissor
Curved

Catalog #	Size
SS-733-16	16 cm (6 1/4")







Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

SS


SCISSORS

Jacobson Micro Scissor
Curved




Catalog #	Size
SS-739-18	18.5 cm (7 1/4")

Adventitia Micro Scissor
Straight, Sharp/Sharp



Catalog #	Size
SS-740-14	14.5 cm (5 3/4")

Micro Scissor
Curved, Sharp/Sharp




Catalog #	Size
SS-741-14	14.5 cm (5 3/4")

Micro Dissecting Scissor
Curved, Blunt/Blunt



Catalog #	Size
SS-743-14	14.5 cm (5 3/4")

Micro Scissor
Straight, Delicate Blades




Catalog #	Size
SS-744-16	16.5 cm (6 1/2")



Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici


Micro Scissor
Curved, Delicate Blades

Catalog #	Size
SS-745-16	16.5 cm (6 1/2")




Micro Scissor
Straight, Extra Delicate Blades

Catalog #	Size
SS-746-16	16.5 cm (6 1/2")




Micro Scissor
Curved, Extra Delicate Blades

Catalog #	Size
SS-747-16	16.5 cm (6 1/2")




Vannas Micro Scissor
Straight

Catalog #	Size
SS-734-16	16 cm (6 1/4")



Vannas Micro Scissor
Curved

Catalog #	Size
SS-735-16	16 cm (6 1/4")





Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

SS

SCISSORS

Millesi Micro Scissor
Straight

Catalog #	Size
SS-736-16	16 cm (6 1/4")




Millesi Micro Scissor
Curved

Catalog #	Size
SS-737-16	16 cm (6 1/4")



Yasargil Micro Scissor
Straight

Catalog #	Size
SS-738-18	18.5 cm (7 1/4")





Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

Micro Coronary Scissor



Catalog #	Size, Profile
SS-751-25	18 cm (7"), Angled 25°
SS-751-45	18 cm (7"), Angled 45°
SS-751-60	18 cm (7"), Angled 60°
SS-751-90	18 cm (7"), Angled 90°
SS-753-00	18 cm (7"), Angled 125°



SS-751-25, Angled 25°



SS-751-45, Angled 45°



SS-751-60, Angled 60°




SS-751-90, Angled 90°



SS-753-00, Angled 125°

Micro Coronary Scissor Titanium



Catalog #	Size, Profile
SS-751-25 T	18 cm (7"), Angled 25°
SS-751-45 T	18 cm (7"), Angled 45°
SS-751-60 T	18 cm (7"), Angled 60°
SS-751-90 T	18 cm (7"), Angled 90°
SS-753-00 T	18 cm (7"), Angled 125°



SS-751-25 T, Angled 25°



SS-751-45 T, Angled 45°



SS-751-60 T, Angled 60°



SS-751-90 T, Angled 90°



SS-753-00 T, Angled 125°



Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici


SS

SCISSORS


Micro Coronary Scissor
Fine Blade



Catalog #	Size, Profile
SS-755-25	16.5 cm (6 1/2"), Angled 25°
SS-755-45	16.5 cm (6 1/2"), Angled 45°
SS-755-60	16.5 cm (6 1/2"), Angled 60°
SS-755-90	16.5 cm (6 1/2"), Angled 90°
SS-757-00	16.5 cm (6 1/2"), Angled 125°




SS-755-25, Angled 25°



SS-755-45, Angled 45°



SS-755-60, Angled 60°



SS-755-90, Angled 90°




SS-757-00, Angled 125°


Micro Coronary Scissor
Fine Blades, One Blade with Probe Tip




Catalog #	Size, Profile
SS-759-25	16.5 cm (6 1/2"), Angled 25°
SS-759-45	16.5 cm (6 1/2"), Angled 45°
SS-759-60	16.5 cm (6 1/2"), Angled 60°
SS-759-90	16.5 cm (6 1/2"), Angled 90°
SS-761-00	16.5 cm (6 1/2"), Angled 125°




SS-759-25, Angled 25°




SS-759-45, Angled 45°



SS-759-60, Angled 60°



SS-759-90, Angled 90°



SS-761-00, Angled 125°




Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

Micro Coronary Scissor
Delicate Blades




Catalog #	Size, Profile
SS-763-25	16.5 cm (6 1/2"), Angled 25°
SS-763-45	16.5 cm (6 1/2"), Angled 45°
SS-763-60	16.5 cm (6 1/2"), Angled 60°
SS-763-90	16.5 cm (6 1/2"), Angled 90°
SS-765-00	16.5 cm (6 1/2"), Angled 125°



Micro Coronary Scissor
Extra Delicate Blades



Catalog #	Size, Profile
SS-767-25	16.5 cm (6 1/2"), Angled 25°
SS-767-45	16.5 cm (6 1/2"), Angled 45°
SS-767-60	16.5 cm (6 1/2"), Angled 60°
SS-767-90	16.5 cm (6 1/2"), Angled 90°
SS-769-00	16.5 cm (6 1/2"), Angled 125°





Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

SS

SCISSORS

Micro Scissor

Straight, Delicate Blades, Round Handle



Catalog #	Size
SS-770-16	16.5 cm (6 1/2")

Micro Scissor

Curved, Delicate Blades, Round Handle



Catalog #	Size
SS-771-16	16.5 cm (6 1/2")

Micro Scissor

Curved, Extra Delicate Blades, Round Handle



Catalog #	Size
SS-773-16	16.5 cm (6 1/2")




Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

Micro Coronary Scissor
Fine Blades, Round Handle



Catalog #	Size, Profile
SS-755-25	16.5 cm (6 1/2"), Angled 25°
SS-755-45	16.5 cm (6 1/2"), Angled 45°
SS-755-60	16.5 cm (6 1/2"), Angled 60°
SS-755-90	16.5 cm (6 1/2"), Angled 90°
SS-777-00	16.5 cm (6 1/2"), Angled 125°



Micro Coronary Scissor
Fine Blades, One Blade with Probe Tip, Round Handle



Catalog #	Size, Profile
SS-779-25	16.5 cm (6 1/2"), Angled 25°
SS-779-45	16.5 cm (6 1/2"), Angled 45°
SS-779-60	16.5 cm (6 1/2"), Angled 60°
SS-779-90	16.5 cm (6 1/2"), Angled 90°
SS-781-00	16.5 cm (6 1/2"), Angled 125°



Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

SS

SCISSORS

Micro Coronary Scissor
Delicate Blades, Round Handle



Catalog #	Size, Profile
SS-783-25	16.5 cm (6 1/2"), Angled 25°
SS-783-45	16.5 cm (6 1/2"), Angled 45°
SS-783-60	16.5 cm (6 1/2"), Angled 60°
SS-783-90	16.5 cm (6 1/2"), Angled 90°
SS-785-00	16.5 cm (6 1/2"), Angled 125°




SS-783-25, Angled 25°



SS-783-45, Angled 45°



SS-783-60, Angled 60°



SS-783-90, Angled 90°



SS-785-00, Angled 125°

Micro Coronary Scissor
Extra Delicate Blades, Round Handle



Catalog #	Size, Profile
SS-787-25	16.5 cm (6 1/2"), Angled 25°
SS-787-45	16.5 cm (6 1/2"), Angled 45°
SS-787-60	16.5 cm (6 1/2"), Angled 60°
SS-787-90	16.5 cm (6 1/2"), Angled 90°
SS-789-00	16.5 cm (6 1/2"), Angled 125°



SS-787-25, Angled 25°



SS-787-45, Angled 45°



SS-787-60, Angled 60°



SS-787-90, Angled 90°




SS-789-00, Angled 125°



Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

Micro Scissor
Straight, Round Handle

Catalog #	Size
SS-790-18	18 cm (7")




Micro Scissor
Curved, Round Handle

Catalog #	Size
SS-791-18	18 cm (7")




Micro Scissor
Curved, One Blade with Probe Tip, Round Handle

Catalog #	Size
SS-793-18	18 cm (7")



Micro Suture Cutting Scissor
Straight, One Toothed Cutting Edge, Round Handle

Catalog #	Size
SS-796-18	18 cm (7")



Micro Scissor
Straight, Round Handle

Catalog #	Size, Blade Size
SS-800-15	15 cm (6"), 10 mm
SS-800-18	18 cm (7"), 10 mm
SS-800-21	21 cm (8 1/4"), 10 mm
SS-800-23	23 cm (9"), 10 mm






Micro Scissor

Straight, Round Handle, Titanium

Catalog #	Size, Blade Size
SS-800-15 T	15 cm (6"), 10 mm
SS-800-18 T	18 cm (7"), 10 mm
SS-800-21 T	21 cm (8 ¼"), 10 mm



Micro Scissor

Curved, Round Handle


Catalog #	Size, Blade Size
SS-801-15	15 cm (6"), 10 mm
SS-801-18	18 cm (7"), 10 mm
SS-801-21	21 cm (8 ¼"), 10 mm
SS-801-23	23 cm (9"), 10 mm



Micro Scissor

Curved, Round Handle, Titanium


Catalog #	Size, Blade Size
SS-801-15 T	15 cm (6"), 10 mm
SS-801-18 T	18 cm (7"), 10 mm
SS-801-21 T	21 cm (8 ¼"), 10 mm



Micro Scissor

Straight, Flat Handle

Catalog #	Size, Blade Size
SS-804-15	15 cm (6"), 10 mm
SS-804-18	18 cm (7"), 10 mm
SS-804-21	21 cm (8 ¼"), 10 mm
SS-804-23	23 cm (9"), 10 mm






Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

Micro Scissor

Straight, Flat Handle, Titanium

Catalog #	Size, Blade Size
SS-804-15 T	15 cm (6"), 10 mm
SS-804-18 T	18 cm (7"), 10 mm
SS-804-21 T	21 cm (8 1/4"), 10 mm



Micro Scissor

Curved, Flat Handle


Catalog #	Size, Blade Size
SS-805-15	15 cm (6"), 10 mm
SS-805-18	18 cm (7"), 10 mm
SS-805-21	21 cm (8 1/4"), 10 mm
SS-805-23	23 cm (9"), 10 mm



Micro Scissor

Curved, Flat Handle, Titanium


Catalog #	Size, Blade Size
SS-805-15 T	15 cm (6"), 10 mm
SS-805-18 T	18 cm (7"), 10 mm
SS-805-21 T	21 cm (8 1/4"), 10 mm



Jacobson Micro Scissor


Straight, Bayonet Shaped

Catalog #	Size
SS-808-23	23 cm (9")






Potts Micro Scissor
Curved, Bayonet Shaped




Catalog #	Size
SS-809-23	23 cm (9")

Micro Scissor
Straight, Bayonet Shaped




Catalog #	Size
SS-810-23	23 cm (9")

Micro Scissor
Curved, Bayonet Shaped




Catalog #	Size
SS-811-23	23 cm (9")

Micro Scissor
Straight, Bayonet Shaped



Catalog #	Size
SS-814-19	18.5 cm (7 1/4")

Micro Scissor
Curved, Bayonet Shaped




Catalog #	Size
SS-815-19	18.5 cm (7 1/4")




Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

Micro Scissor
Straight, Bayonet Shaped




Catalog #	Size
SS-816-19	19 cm (7 1/2")
SS-816-21	21.5 cm (8 1/2")

Micro Scissor
Curved Downwards, Bayonet Shaped



Catalog #	Size
SS-817-19	18.5 cm (7 1/4")
SS-817-21	21.5 cm (8 1/2")

Micro Scissor
Curved Upwards, Bayonet Shaped




Catalog #	Size
SS-819-19	18.5 cm (7 1/4")
SS-819-21	21.5 cm (8 1/2")

Jaboma Micro Scissor
Angled, One Blade Probe Tip



Catalog #	Size, Blade Size
SS-825-18	18 cm (7"), 12 mm

Jaboma Micro Scissor
Angled, One Blade Probe Tip



Catalog #	Size, Blade Size
SS-827-18	18 cm (7"), 12 mm



SURTEX
INSTRUMENTS

Crafted to Excellence

Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici

SS

SCISSORS

Jaboma Micro Scissor (For Cleft Palate)
Angled

Catalog #	Size, Blade Size
SS-829-18	18 cm (7"), 12 mm






Operating Scissors
Chirurgische Scheren
Tijeras quirúrgicas
Ciseaux chirurgicaux
Forbici per chirurgia

UltraCut™ TC Operating Scissor
Straight, Blunt/Blunt

ULTRA.Cut™



Catalog #	Size
SS-850-14	14.5 cm (5 ¾")



UltraCut™ TC Operating Scissor
Curved, Blunt/Blunt

ULTRA.Cut™


Catalog #	Size
SS-851-14	14.5 cm (5 ¾")



UltraCut™ TC Operating Scissor
Straight, Sharp/Blunt

ULTRA.Cut™


Catalog #	Size
SS-852-14	14.5 cm (5 ¾")



XTSCut™ TC Operating Scissor
Straight, Sharp/Blunt

XTS.Cut™



Catalog #	Size
SS-852-14 SC	14.5 cm (5 ¾")



UltraCut™ TC Operating Scissor
Curved, Sharp/Blunt

ULTRA.Cut™

Catalog #	Size
SS-853-14	14.5 cm (5 ¾")





Operating Scissors
Chirurgische Scheren
Tijeras quirúrgicas
Ciseaux chirurgicaux
Forbici per chirurgia


SS

SCISSORS

XTSCut™ TC Operating Scissor
Curved, Sharp/Blunt




Catalog #	Size
SS-853-14 SC	14.5 cm (5 ¾")



UltraCut™ TC Operating Scissor
Straight, Sharp/Sharp



Catalog #	Size
SS-854-14	14.5 cm (5 ¾")



UltraCut™ TC Operating Scissor
Curved, Sharp/Sharp



Catalog #	Size
SS-855-14	14.5 cm (5 ¾")



UltraCut™ TC Mayo Dissecting Scissor
Straight



Catalog #	Size
SS-860-14	14.5 cm (5 ¾")
SS-860-17	17 cm (6 ¾")
SS-860-23	23 cm (9")



XTSCut™ TC Mayo Dissecting Scissor
Straight



Catalog #	Size
SS-860-14 SC	14.5 cm (5 ¾")
SS-860-17 SC	17 cm (6 ¾")
SS-860-23 SC	23 cm (9")






Operating Scissors
Chirurgische Scheren
Tijeras quirúrgicas
Ciseaux chirurgicaux
Forbici per chirurgia

UltraCut™ TC Mayo Dissecting Scissor
Curved

ULTRA Cut™

Catalog #	Size
SS-861-14	14.5 cm (5 ¾")
SS-861-17	17 cm (6 ¾")
SS-861-23	23 cm (9")



XTSCut™ TC Mayo Dissecting Scissor
Curved

XTS Cut™

Catalog #	Size
SS-861-14 SC	14.5 cm (5 ¾")
SS-861-17 SC	17 cm (6 ¾")
SS-861-23 SC	23 cm (9")



UltraCut™ TC Mayo-Stille Dissecting Scissor
Straight

ULTRA Cut™

Catalog #	Size
SS-890-15	15 cm (6")
SS-890-17	17 cm (6 ¾")



XTSCut™ TC Mayo-Stille Dissecting Scissor
Straight

XTS Cut™

Catalog #	Size
SS-890-15 SC	15 cm (6")
SS-890-17 SC	17 cm (6 ¾")





Operating Scissors
Chirurgische Scheren
Tijeras quirúrgicas
Ciseaux chirurgicaux
Forbici per chirurgia


SS

SCISSORS

UltraCut™ TC Mayo-Stille Dissecting Scissor
Curved

ULTRACut™


Catalog #	Size
SS-891-15	15 cm (6")
SS-891-17	17 cm (6 ¾")



XTSCut™ TC Mayo-Stille Dissecting Scissor
Curved

XTSCut™

Catalog #	Size
SS-891-15 SC	15 cm (6")
SS-891-17 SC	17 cm (6 ¾")



UltraCut™ TC Metzenbaum Dissecting Scissor
Straight

ULTRACut™

Catalog #	Size
SS-900-11	11.5 cm (4 ½")



UltraCut™ TC Metzenbaum Dissecting Scissor
Curved

ULTRACut™

Catalog #	Size
SS-901-11	11.5 cm (4 ½")



UltraCut™ TC Metzenbaum Dissecting Scissor
Straight

ULTRACut™

Catalog #	Size
SS-902-14	14.5 cm (5 ¾")





Operating Scissors
Chirurgische Scheren
Tijeras quirúrgicas
Ciseaux chirurgicaux
Forbici per chirurgia

XTSCut™ TC Metzenbaum Dissecting Scissor
Straight



Catalog #	Size
SS-902-14 SC	14.5 cm (5 ¾")



UltraCut™ TC Metzenbaum Dissecting Scissor
Curved




Catalog #	Size
SS-903-14	14.5 cm (5 ¾")



XTSCut™ TC Metzenbaum Dissecting Scissor
Curved



Catalog #	Size
SS-903-14 SC	14.5 cm (5 ¾")



UltraCut™ TC Metzenbaum Dissecting Scissor
Straight



Catalog #	Size
SS-904-18	18 cm (7")
SS-904-20	20.5 cm (8")
SS-904-23	23 cm (9")
SS-904-25	26 cm (10 ¼")
SS-904-28	28.5 cm (11 ¼")
SS-904-30	31 cm (12 ¼")





Operating Scissors
Chirurgische Scheren
Tijeras quirúrgicas
Ciseaux chirurgicaux
Forbici per chirurgia

SS



SCISSORS

XTSCut™ TC Metzenbaum Dissecting Scissor
Straight




Catalog #	Size
SS-904-18 SC	18 cm (7")
SS-904-20 SC	20.5 cm (8")
SS-904-23 SC	23 cm (9")

UltraCut™ TC Metzenbaum Dissecting Scissor
Curved



Catalog #	Size
SS-905-18	18 cm (7")
SS-905-20	20.5 cm (8")
SS-905-23	23 cm (9")
SS-905-25	26 cm (10 1/4")
SS-905-28	28.5 cm (11 1/4")
SS-905-30	31 cm (12 1/4")

XTSCut™ TC Metzenbaum Dissecting Scissor
Curved



Catalog #	Size
SS-905-18 SC	18 cm (7")
SS-905-20 SC	20.5 cm (8")
SS-905-23 SC	23 cm (9")

UltraCut™ TC Metzenbaum -Thorek Dissecting Scissor
Curved



Catalog #	Size
SS-907-20	20 cm (8")
SS-907-23	23 cm (9")
SS-907-28	28 cm (11")



Operating Scissors, Dissecting Scissors
 Chirurgische Scheren, Präparierscheren
 Tijeras quirúrgicas, Tijeras de disección
 Ciseaux chirurgicaux, Ciseaux à dissection
 Forbici per chirurgia, Forbici per dissezione

UltraCut™ TC Joseph Dissecting Scissor
 Straight

ULTRACut™


Catalog #	Size
SS-910-14	14.5 cm (5 3/4")



UltraCut™ TC Joseph Dissecting Scissor
 Curved

ULTRACut™

Catalog #	Size
SS-911-14	14.5 cm (5 3/4")



UltraCut™ TC Lexer Dissecting Scissor
 Straight

ULTRACut™

Catalog #	Size
SS-914-16	16 cm (6 1/4")
SS-914-21	21 cm (8 1/4")



XTSCut™ TC Lexer Dissecting Scissor
 Straight

XTSCut™


Catalog #	Size
SS-914-16 SC	16 cm (6 1/4")



UltraCut™ TC Lexer Dissecting Scissor
 Curved

ULTRACut™

Catalog #	Size
SS-915-16	16 cm (6 1/4")
SS-915-21	21 cm (8 1/4")





Operating Scissors, Dissecting Scissors
Chirurgische Scheren, Präparierscheren
Tijeras quirúrgicas, Tijeras de disección
Ciseaux chirurgicaux, Ciseaux à dissection
Forbici per chirurgia, Forbici per dissezione

XTSCut™ TC Lexer Dissecting Scissor
Curved



Catalog #	Size
SS-915-16 SC	16 cm (6 1/4")



UltraCut™ TC Mayo-Lexer Dissecting Scissor
Curved



Catalog #	Size
SS-917-16	16 cm (6 1/4")



UltraCut™ TC Metzenbaum-Fine Delicate Dissecting Scissor
Straight, Slender Pattern




Catalog #	Size
SS-920-14	14.5 cm (5 3/4")



XTSCut™ TC Metzenbaum-Fine Delicate Dissecting Scissor
Straight, Slender Pattern




Catalog #	Size
SS-920-14 SC	14.5 cm (5 3/4")



UltraCut™ TC Metzenbaum-Fine Delicate Dissecting Scissor
Curved, Slender Pattern



Catalog #	Size
SS-921-14	14.5 cm (5 3/4")





Delicate Dissecting Scissors
Feine Präparierscheren
Tijeras finas de disección
Ciseaux fins à dissection
Forbici fini per dissezione

XTSCut™ TC Metzenbaum-Fine Delicate Dissecting Scissor
Curved, Slender Pattern




Catalog #	Size
SS-921-14 SC	14.5 cm (5 ¾")



UltraCut™ TC Metzenbaum-Fine Delicate Dissecting Scissor
Straight, Slender Pattern




Catalog #	Size
SS-922-18	18 cm (7")
SS-922-20	20.5 cm (8")
SS-922-23	23 cm (9")
SS-922-25	25 cm (9 ¾")
SS-922-28	28 cm (11")



XTSCut™ TC Metzenbaum-Fine Delicate Dissecting Scissor
Straight, Slender Pattern




Catalog #	Size
SS-922-18 SC	18 cm (7")
SS-922-20 SC	20.5 cm (8")
SS-922-23 SC	23 cm (9")



UltraCut™ TC Metzenbaum-Fine Delicate Dissecting Scissor
Curved, Slender Pattern



Catalog #	Size
SS-923-18	18 cm (7")
SS-923-20	20.5 cm (8")
SS-923-23	23 cm (9")
SS-923-25	25 cm (9 ¾")
SS-923-28	28.5 cm (11 ¼")





Delicate Dissecting Scissors
Feine Präparierscheren
Tijeras finas de disección
Ciseaux fins à dissection
Forbici fini per dissezione

SS


SCISSORS

XTSCut™ TC Metzenbaum-Fine Delicate Dissecting Scissor
Curved, Slender Pattern




Catalog #	Size
SS-923-18 SC	18 cm (7")
SS-923-20 SC	20.5 cm (8")
SS-923-23 SC	23 cm (9")

UltraCut™ TC Metzenbaum-Fine Delicate Dissecting Scissor
Straight, Sharp




Catalog #	Size
SS-926-14	14.5 cm (5 3/4")

UltraCut™ TC Metzenbaum-Fine Delicate Dissecting Scissor
Curved, Sharp



Catalog #	Size
SS-927-14	14.5 cm (5 3/4")

UltraCut™ TC Metzenbaum-Fine Delicate Dissecting Scissor
Straight, Sharp



Catalog #	Size
SS-928-18	18 cm (7")
SS-928-20	20 cm (8")

UltraCut™ TC Metzenbaum-Lahey Dissecting Scissor
Curved, Slender Pattern



Catalog #	Size
SS-931-14	14.5 cm (5 3/4")




Delicate Dissecting Scissors, Nasal Scissors
 Feine Präparierscheren, Nasenscheren
 Tijeras finas de disección, Tijeras nasales
 Ciseaux fins à dissection, Ciseaux nasaux
 Forbici fini per dissezione, Forbici nasali

UltraCut™ TC Toennis-Adson Dissecting Scissor
 Curved, Slender Pattern

ULTRACut™


Catalog #	Size
SS-933-17	17.5 cm (7")



UltraCut™ TC Metzenbaum-Lahey Dissecting Scissor
 Curved, Slender Pattern

ULTRACut™


Catalog #	Size
SS-935-20	20.5 cm (8")
SS-935-23	23 cm (9")



UltraCut™ TC Metzenbaum-Fine Delicate Dissecting Scissor
 Curved, S Shaped, Slender Pattern

ULTRACut™

Catalog #	Size
SS-937-18	18 cm (7")
SS-937-20	20.5 cm (8")
SS-937-23	23 cm (9")



UltraCut™ TC Fomon Nasal Scissor

ULTRACut™

Catalog #	Size
SS-939-15	14.5 cm (5 ¾")



UltraCut™ TC Cottle Nasal Scissor

ULTRACut™

Catalog #	Size
SS-941-16	16 cm (6 ¼")



Face-lift Scissors
Face-lift-Scheren
Tijeras para Face-Lift
Ciseaux pour Face-lift
Forbici per Face-lift

SS
SCISSORS

UltraCut™ TC Gorney Face-Lift Scissor
Straight, One Toothed Cutting Edge

ULTRACut™

Catalog #	Size
SS-944-20	20 cm (8")
SS-944-23	23 cm (9")



UltraCut™ TC Rees Face-Lift Scissor
Curved, One Toothed Cutting Edge

ULTRACut™


Catalog #	Size
SS-947-17	17 cm (6 3/4")



UltraCut™ TC Gregory Face-Lift Scissor
One Toothed Cutting Edge

ULTRACut™

Catalog #	Size
SS-949-14	14 cm (5 1/2")



UltraCut™ TC Kaye Face-Lift Scissor
One Toothed Cutting Edge

ULTRACut™


Catalog #	Size
SS-951-15	15 cm (6")
SS-951-18	18 cm (7")



UltraCut™ TC Stella-S Face-Lift Scissor
S Shaped, One Toothed Cutting Edge

ULTRACut™

Catalog #	Size
SS-953-18	18 cm (7")







Face-lift Scissors
 Face-lift-Scheren
 Tijeras para Face-Lift
 Ciseaux pour Face-lift
 Forbici per Face-lift

UltraCut™ TC Gorney Face-Lift Scissor
 One Toothed Cutting Edge

ULTRACut™



Catalog #	Size
SS-959-18	18 cm (7")



UltraCut™ TC Freeman-Kaye Face-Lift Scissor
 One Toothed Cutting Edge

ULTRACut™


Catalog #	Size
SS-961-18	18 cm (7")
SS-963-23	23 cm (9")



UltraCut™ TC Wertheim Gynecological Scissor
 Curved

ULTRACut™

Catalog #	Size
SS-965-23	23 cm (9")
SS-965-35	35 cm (14")



UltraCut™ TC Sims Gynecological Scissor
 Straight

ULTRACut™

Catalog #	Size
SS-968-23	23 cm (9")



UltraCut™ TC Sims Uterine Scissor
 Curved

ULTRACut™

Catalog #	Size
SS-969-23	23 cm (9")





Hysterectomy Scissors
Hysterektomiescheren
Tijeras para histerectomía
Ciseaux à hystérectomie
Forbici per isterectomia

SS

SCISSORS

UltraCut™ TC Parametrium Hysterectomy Scissor
Strongly Curved

ULTRACut™


Catalog #	Size
SS-971-01	22.5 cm (8 ¾")



UltraCut™ TC Parametrium Hysterectomy Scissor
Curved

ULTRACut™


Catalog #	Size
SS-971-02	23 cm (9")



UltraCut™ TC Stevens Tenotomy Scissor
Straight, Blunt/Blunt

ULTRACut™


Catalog #	Size
SS-972-10	10.5 cm (4 ¼")
SS-972-11	11.5 cm (4 ½")



UltraCut™ TC Stevens Tenotomy Scissor
Curved, Blunt/Blunt

ULTRACut™


Catalog #	Size
SS-973-10	10.5 cm (4 ¼")
SS-973-11	11.5 cm (4 ½")



UltraCut™ TC Stevens Tenotomy Scissor
Straight, Sharp/Sharp

ULTRACut™

Catalog #	Size
SS-974-10	10.5 cm (4 ¼")
SS-974-11	11.5 cm (4 ½")






Iris Scissors
Irissscheren
Tijeras para iridectomía
Ciseaux à iridectomie
Forbici per iridectomia

UltraCut™ TC Stevens Tenotomy Scissor
Curved, Sharp/Sharp

ULTRA Cut™


Catalog #	Size
SS-975-10	10.5 cm (4 1/4")
SS-975-11	11.5 cm (4 1/2")



UltraCut™ TC Goldman-Fox Gum Scissor
Straight

ULTRA Cut™


Catalog #	Size
SS-976-13	13 cm (5")



UltraCut™ TC Goldman-Fox Gum Scissor
Curved

ULTRA Cut™

Catalog #	Size
SS-977-13	13 cm (5")



UltraCut™ TC Iris Scissor
Straight

ULTRA Cut™

Catalog #	Size
SS-978-11	11.5 cm (4 1/2")



XTSCut™ TC Iris Scissor
Straight

XTS Cut™

Catalog #	Size
SS-978-11 SC	11.5 cm (4 1/2")





Iris Scissors
Irisscheren
Tijeras para iridectomía
Ciseaux à iridectomie
Forbici per iridectomia

SS

SCISSORS

UltraCut™ TC Iris Scissor
Curved

ULTRA Cut™

Catalog #	Size
SS-979-11	11.5 cm (4 ½")



XTSCut™ TC Iris Scissor
Curved

XTS Cut™

Catalog #	Size
SS-979-11 SC	11.5 cm (4 ½")




UltraCut™ TC DeBakey Vascular Scissor


ULTRA Cut™

Catalog #	Size
SS-981-25	19 cm (7 ½")
SS-981-45	18.5 cm (7 ¼")
SS-981-60	18 cm (7")


Catalog #	Size
SS-983-25	23 cm (9")
SS-983-45	23 cm (9")
SS-983-60	22 cm (8 ¾")



SS-981-25, Angled 25°




SS-981-45, Angled 45°




SS-981-60, Angled 60°



SS-983-25, Angled 25°



SS-983-45, Angled 45°



SS-983-60, Angled 60°




Iris Scissors
 Irisscheren
 Tijeras para iridectomía
 Ciseaux à iridectomie
 Forbici per iridectomia

UltraCut™ TC Iris Scissor
 Angled

ULTRA.Cut™



Catalog #	Size
SS-985-11	11 cm (4 1/2")



UltraCut™ TC Ragnell (Kilner) Dissecting Scissor
 Curved

ULTRA.Cut™



Catalog #	Size
SS-987-12	12.5 cm (5")
SS-987-15	15 cm (6")
SS-987-18	18 cm (7")



UltraCut™ TC Kelly Dissecting Scissor
 Straight

ULTRA.Cut™



Catalog #	Size
SS-988-16	16 cm (6 1/4")



UltraCut™ TC Kelly Dissecting Scissor
 Curved

ULTRA.Cut™



Catalog #	Size
SS-989-16	16 cm (6 1/4")



UltraCut™ TC Universal Wire Cutting Scissor
 Angled One Toothed Cutting Edge

ULTRA.Cut™

Catalog #	Size
SS-991-12	12.5 cm (5")





Dissecting Scissors, Gum Scissors
Präparierscheren, Zahnfleischscheren
Tijeras de disección, Tijeras para encías
Ciseaux à dissection, Ciseaux à gencives
Forbici per dissezione, Forbici per gengiva


SS

SCISSORS

UltraCut™ TC Dean Tonsil Scissor
S Shaped, One Toothed Cutting Edge

ULTRACut™

Catalog #	Size
SS-993-17	17.5 cm (7")



UltraCut™ TC Locklin Gum Scissor
Angled, One Toothed Cutting Edge

ULTRACut™

Catalog #	Size
SS-995-16	16 cm (6 1/4")



UltraCut™ TC Locklin Gum Scissor
Angled, S Shaped, One Toothed Cutting Edge

ULTRACut™

Catalog #	Size
SS-997-16	16 cm (6 1/4")



UltraCut™ TC Lister Bandage Scissor
Angled

ULTRACut™

Catalog #	Size
SS-999-18	18 cm (7")



